
101 277

समुदायमा आधारित पानीको गुणस्तर अनुगमन तथा
व्यवस्थापन पुस्तिका
**Community Based Water Quality Monitoring and
Drinking Water Management Manual**

Edited by:

Dr. Durga Datt Joshi
Mr. Andres Sanchez
Mr. Banrey Dutka
Dr. Madhav Karki
Mr. Mahendra Maharjan

Published by:

National Zoonoses and Food Hygiene Research Centre (NZFHRC)

Supported by:

International Development Research Centre (IDRC), Ottawa, Canada

2004

List of Publications of NZFHRC with ISBN Number

1. Present Situation on the Control of Human and Canine Rabies In Nepal. ISBN No. 99933 - 29 - 35 - 7
2. Make Milk Not Methane Gas In Nepal. ISBN No. 99933-29-33-9.
3. Water Quality Status of Kathmandu Aquatox 2000 Programme in Nepal. ISBN No. 99933-29-34-7.
4. Food Safety in Nepal. ISBN No. 99933-29-32-0
5. Cost Analysis of Milk Production Cow and Buffalo in Nepal. ISBN No. 99933-29-27-4
6. Economic Impact Evaluation of Milk Producers in Mid Western Milk Supply Scheme of Dair Development Corporation. ISBN No. 99933-29-25-8
7. Epidemiology of Cystic Echinococcosis/Hydatidosis Distribution and Transmission Patterns i Kathmandu Nepal. ISBN No. 99933-29- 21-5
8. Public Health and Medical Officers Training Snakebite Management in Nepal. ISBN No. 99933-29- 23-1
9. Epidemiological Surveillance of Human and animal Brucellosis in Milk Collection Area of Dairy Development Corporation in Nepal. ISBN No. 99933-29-20-7
10. Yak and Chauri Husbandry in Nepal. ISBN No. 99933-29-09-6
11. Urban Echinococcosis in Health Transition in Nepal. ISBN No. 99933-29-08-8
12. Marketing Status of Milk and Milk Products in Kathmandu, Lalitpur and Bhaktapur. ISBN No. 9993- 29-11-8
13. Inter-Regional Workshop on Anthrax Proceedings. ISBN No. 99933-29-06-1
14. Rabies Present Status and Future Strategy to Control in Animals and Human, Nepal. ISBN No. 99933 29-05-3.
15. Yak (Bos Grunniens) in Central Asia. ISBN No. 99933-29-03-7
16. Animal Slaughtering and Meat Marketing Practices in Nepal. ISBN No. 99933-29-02-9.
17. An Assessment of the Yak Cheese Factories in Nepal. ISBN No. 99933-29- 00-2
18. Urban Ecosystem Health Status in Ward 19 and 20 of Kathmandu. ISBN No. 99933-29-01-0
19. Animal Health, Livestock and Dairy Development of Nepal. ISBN No. 99933-546-1-2
20. Livestock Trade Statistics of Nepal. . ISBN No. 99933-546-2-7
21. Veterinary Public Health Hazards in Nepal. ISBN No. 99933-546-0-0
22. Report on Epidemiology of Rabies in Nepal Year 1998-1999. ISBN No. 99933-29-16-9
23. Epidemiological Surveillance Report on Snakebite in Nepal for the Years 2055 to 2056 B.S. ISBN No. 99933-29-17-7
24. Epidemiological Surveillance Report on Japanese Encephalitis in Nepal. ISBN No. 99933-29-1-3
25. Epidemiological Surveillance Report on Visceral Leishmaniasis for the Year 2055/2056 B.S. ISBN No 99933-29-18-5
26. Humane Slaughtering Management Training W-S. ISBN No. 99933-49-40-2
27. Humane Slaughtering Management Training. ISBN No. 99933-49-41-0
28. Community Based Water Quality Monitoring and Drinking Water Management Training for Selected Ward Members of KMC (Nepali + English). ISBN No. 99933-49-42-9
29. Community Based Water Quality Monitoring and Drinking Water Management Sensitization Workshop and Training A Pilot Project in Wards 19 and 20 of KMC Nepal. ISBN No. 99933-49-43-7
30. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Kirtipur Municipality, Kathmandu District. ISBN No. 99933-49-44-5
31. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Kathmandu Metropolitan City, Kathmandu District. ISBN No. 99933-49-45-3
32. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Lalitpur Sub-Metropolitan City, Lalitpur District. ISBN No. 99933-49-46-1
33. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Madhyapur Thimi Municipality, Bhaktapur District. ISBN No. 99933-49-47-X
34. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Bhaktapur Municipality, Bhaktapur District. ISBN No. 99933-49-48-8
35. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Banepa Municipality, Kavre District. ISBN No. 99933-49-49-6
36. Final Report on Dog Rabies Vaccination and Mass Awareness Campaign in Panauti Municipality, Kavre District. ISBN No. 99933-49-50-X

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**Community Based Water Quality Monitoring and
Drinking Water Management Manual**

Dr. Durga Datt Joshi

Director

National Zoonoses & Food Hygiene Research Centre

Mr. Andres Sanchez

Programme Officer

International Development Research Centre (IDRC), Ottawa, Canada

Mr. Barney Dutka

Water Consultant

International Development Research Centre (IDRC), Ottawa, Canada

Dr. Madhav Karki

Regional Advisor

International Development Research Centre (IDRC), New Delhi, India.

Mr. Mahendra Maharjan

Microbiologist

National Zoonoses & Food Hygiene Research Centre

Published by:

National Zoonoses and Food Hygiene Research Centre

Supported by:

International Development Research Centre (IDRC), Ottawa, Canada.

2004

663.6 01/12/04

PUBLISHER:

National Zoonoses and Food Hygiene Research Centre (NZFHRC)

P.O.Box: 1885, Tahachal, Kathmandu, Nepal

Telephone: +977-1-4270667, Fax: +977-1-4272694

E-mail: ddjoshi@healthnet.org.np

Website: www.nzfhr.org.np

ISBN No.: 99933 - 29 - 38 - X

First Edition Publication in 2001.

Second Edition Publication in 2004.

© Copyright NZFHRC 2001/04.

A Catalogue record for this book is available from the Central Library Tribhuvan University Kirtipur Kathmandu Nepal.

No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, and recording or otherwise, without the permission of the publisher. The designations employed and the presentation of the material in this publication do not imply the expression of any opinion whatsoever on the part of the Secretariat of the NZFHRC.

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TABLE OF CONTENTS

Conversion Table from Nepali Calendar Year to Gregorian Calendar year. B.S to A.D

Preface and Acknowledgement

1.	Organization of Manual	2
1.1	Introduction	2
1.2	Notes to the Facilitator	6
1.3	Proposed Training Schedules	8
2.	Collecting Basic Information	10
2.1	Community Mapping	10
2.2	Water Use Survey	14
2.3	Transect Walk	18
3.	Monitoring Bacterial Water Quality	22
3.1	How Water Becomes Contaminated	24
3.2	Bacterial Indicator Systems	32
3.3	The H ₂ S Test	38
3.4	Sample Collection & Processing	50
3.5	Coliplate & Colistrip Tests	56
3.6	Recording and Interpreting Results	58
3.7	Cleaning & Disposing of contaminated samples	64
4.	Setting Up A Local Water Quality Monitoring Program	66
4.1	Identifying Sampling Points	68
4.2	Roles and Responsibilities	76
4.3	Planning the Program	84
4.4	Reporting Requirements	88
5.	Some Remedial Actions	94
5.1	Filtration	96
5.2	Disinfection	99
5.3	Storage	115
6.	Training Evaluation & Closing	120
Annexes:	Trainee Handouts	124
	Community Mapping	125
	Water Use Survey	126
	Transect Walk	128
	The H ₂ S Test	134
	Coliplate™ & Colistrip™ Test Guide	151

Module 7: Water Quality Control Techniques (IDRC)

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Conversion Table from Nepali Calendar Year to Gregorian Calendar year. B.S to A.D

Nepali Month Conversion Table	
Nepali Calendar Month	Gregorian Calendar Equivalent
Baishak	Middle of April to middle of June
Jestha	Middle of May to middle of July
Asadh	Middle of June to middle of July
Srawan	Middle of July to middle of August
Bhadra	Middle of August to middle of September
Aswin	Middle of September to middle of October
Kartik	Middle of October to middle of November
Mangsir	Middle of November to middle of December
Poush	Middle December of to middle of January
Magh	Middle of January to middle of February
Falgun	Middle of February to middle of March
Chaitra	Middle of March to middle of April

Nepali Year Conversion Table	
Nepali Calendar Year (Begins in Mid April) B.S. (Bikram Samvat)	Gregorian Calendar Year A.D.
2045	1988/89
2046	1989/90
2047	1990/91
2048	1991/92
2049	1992/93
2050	1993/94
2051	1994/95
2052	1995/96
2053	1996/97
2054	1997/98
2055	1998/99
2056	1999/2000
2057	2000/2001
2058	2001/2002
2059	2002/2003
2060	2003/2004

ACKNOWLEDGMENTS

The materials presented in this training guide have been compiled and adapted from a number of different sources listed below. Their contribution is gratefully acknowledged.

Donald Sharp (1997). Training Program: Water Quality Monitoring. Document prepared for the Water Supply and Sanitation for Low Income Communities Project (Indonesia).

Bernard Dutka (1996). Bacteriological Methods Manual for Testing Drinking Water in Rural and Isolated Indonesian Villages. Environment Canada. National Water Research Institute. Burlington, Ontario, Canada.

Bernard Dutka and Andrés Sánchez (1998). Module 7: Water Quality Control Techniques. Lorena Aguilar Revelo, Elías Rosales Escalante, Editors. First Edition. Fundación Tecnológica de Costa Rica. Translated into English from Spanish original.

Lorena Aguilar Revelo, Elías Rosales Escalante, Editors (1998). Module 5: Slow Sand Filters. First Edition. Fundación Tecnológica de Costa Rica. Translated into English from Spanish original.

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1. Organization of Manual

1.1 Introduction

Through out the developing world, limited government resources, aging infrastructure and rapid urbanization lead to stressed water treatment and supply systems that fail to serve the basic needs of the poor. Cities typically experience rapid increases in disadvantaged populations that settle in informal or unplanned areas with inadequate water and sanitation services. In these situations, the responsibility for making drinking water safe for human consumption often falls by default to community residents.

The community-based Water Quality Monitoring (WQM) approach promoted through these training aims to develop and test an appropriate and affordable monitoring alternative for the urban and rural government that promotes collaboration between village/city community organizations, NGOs and water and health authorities. Its goal is to increase monitoring coverage and lead to a better understanding by the population of the interrelationships between water, health and prevention of diarrheal diseases. The approach will involve community groups in planning the monitoring their own water supply and sanitation systems, collect water samples, perform simple water quality tests and make decisions on remedial actions in collaboration with local and city authorities. In order to establish such a program, communities must first understand the need for WQM, be willing to contribute to the success of the program and be prepared to take action to correct problems when identified.

To insure that a WQM program is sustainable, community members must be involved in its design from the start. These training materials and water quality tests are designed for use in local wards. They are intended to take full advantage of participatory approaches to obtain the information needed for establishing the program while at the same time develop skills among those who will be implementing it:

The objectives of this training program will be as follows:

1. To provide information and hands-on training that will motivate community people, organizations, NGOs and water and health authorities to establish a WQM program.
2. To determine the environmental health conditions that exist in each ward enrolled in the program, including perceptions of water quality, and water use practices.
3. To develop the technical skills and knowledge needed to carry out the program.
4. To define roles and responsibilities, and design an implementation plan suitable to addressing the needs of the participating wards.

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The program will impart to community members and local organizations **new knowledge** and develop the **necessary skills** needed to establish a ward-level WQM program; and, will **train trainers** in the use of techniques intended for replicating the program in other communities of the city and village.

The first step in implementing this program is to test and adapt the training materials contained herein. The topics covered may be grouped into four categories.

- “Community Mapping”, “Water Use Survey”, and “Transect Walk” sections involve the collection of background knowledge and are meant to establish the need for a WQM program in a particular ward.
- The section “Monitoring Bacterial Water Quality” is technical in nature and is meant for those who will be implementing the program.
- “Setting Up A Local Water Quality Monitoring Program” covers organizational and administrative aspects and is intended for those who will be planning and managing the program. It also addresses policy issues that cannot be finalized until after consultation with all concerned partners.
- Finally, section 5: “Some Remedial Actions” presents simple alternatives for households and communities that lack potable water, allowing them to protect themselves against a variety of waterborne pathogens and reduce the incidence of waterborne diarrheal diseases.

As these modules involve Hygiene and Sanitation Education (HSE), it is recommended that where appropriate, the content be fully integrated into the ongoing community HSE programs.

It is also important to note that the approaches and techniques being introduced here have been simplified for local-level use. Wherever possible an attempt has been made to reduce the number of steps involved in performing a task, and to use locally available resources.

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1.2 Notes to the Facilitator

The first step in this exercise is to test the training materials presented here. This will be done through the technology transfer training workshops and pilot ward projects. It is suggested that not more than 20 - 25 persons take part in any one training exercise at any one time. Large numbers make it difficult for all to participate and are harder to supervise.

The teaching techniques used in this program are meant to involve the participants in the learning process as much as possible: i.e. the "hands-on" approach. This means that it is important that the answers to key questions come as much as possible from the participants not the facilitator. In other words, it is anticipated that learning will take place through a process of discovery and active participation, not through lecture and passive observation. For this method to be effective it is necessary that the facilitator motivate the participants to take part in the exercises. This is done by making the presentation as interesting as possible. Telling jokes, songs and role playing help to relax people so that they don't feel shy. Films, posters and handout materials can help to reinforce the messages being introduced.

At the start of each day's session, the activities of the previous day are reviewed by asking a participant to summarize what was done and what was learned. Each day the participants are reminded about the purpose of the program so that they do not lose sight of the overall objective of the training. If a section or procedure was not fully understood, these can be reviewed and/or repeated. Also, the facilitator should clearly show the linkage between each section during the review period. In this way, important messages are continually being reinforced. At the close of each day's session, time is reserved for questions, answers and comments.

It should be emphasized in the beginning that there are no right or wrong answers and that all questions are welcome. In fact, when asking questions, if a wrong answer is given, the facilitator does not point out the error. Instead the facilitator says nothing but continues asking the question until the correct answer is given. He/she then reinforces the correct answer by getting the participant to repeat it for everyone to hear. If the majority of the audience is literate the answers can be summarized on a flip chart or chalkboard for everyone to see. Pictures (sketches and photos) also help. When showing what a good situation looks like, a set of photographs can be used. Ideally, someone in the group, or an artist hired by the project, can sketch certain situations relevant to the community as no two communities are identical.

The venue is also important. As it may be necessary to make special arrangements for a room or place that can accommodate all of the participants, it is important that preparations and visits to the community be made well in advance of the scheduled meeting to avoid disappointing results. The venue, whether inside or outside should be at a place where there are few distractions, where it is well lit and where there is ample space for flip charts and for demonstrations. Also, as much of the training is on-the-job type training, it is important that sufficient time is allocated to finding suitable nearby examples for field work.

Finally, no one likes to come to a training exercise without refreshments. Provision for food and drink should be arranged in advance.

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1.3 Proposed Training Schedules

Time	Tech. Transfer Training
Day 1: AM	Introductions
PM	Workshop objectives Overview: collecting basic information 3.1 How water becomes contaminated 3.2 Bacterial Indicator Systems
Day 2: AM	3.3 The H ₂ S test (Lab)
PM	Preparation of paper strips for ward pilot programs
Day 3: AM	3.4 Sample collection & processing
PM	3.5 Coliplate & Colistrip tests (Lab) Overview: Pilot WQM Programs
Day 4: AM	3.6 Recording & interpreting results
PM	3.7 Cleaning & disposing of contaminated samples Overview: Remedial Actions Training Evaluation & Closing
Day 5: AM	
PM	Preparing Ward Kits

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2. Collecting Basic Information

2.1 Community Mapping

- Purpose:**
- 1) To document the layout of the ward, identify important landmarks and locate the position of the existing water supply and sanitation facilities.
 - 2) To promote an awareness among the villagers of their existing environmental health conditions, both good and bad.

Proposed Participants: The community at large (men, women and children), the water user's committee (if any), relevant NGO and government personnel.

Expected Outcome: A ward map will be produced showing the location of all existing water supply and sanitation facilities (latrines, wells, springs, pipelines, water sources, etc.) in relationship to key landmarks. This map will later be used to plot the location of potential sources of pollution and determine water sampling points.

Materials:

- ◆ Large sheets of paper at least 1m square (eg. flip chart paper)
- ◆ Masking tape
- ◆ Coloured pens, markers or crayons
- ◆ Pads, pens and rulers for the participants to take notes

Timing: 2 x 45 min.

Methods: (should be adapted to the experience/ preference of the facilitator and community context)

1. After a brief explanation of the purpose of the training program and an introduction to those present, the facilitator asks for a volunteer to come to the front where the paper has been put on the wall, to make a map of the community (village, ward and city).
2. The facilitator asks the volunteer to choose a coloured marker and draw in the main roads and any streams, rivers canals or ponds that run through the community. The audience can help.
3. Next, the facilitator asks another volunteer to draw in the religious centres of worship, health clinics, schools, and other community centres.
4. If the community is small, the location of other main buildings and public squares can be included.

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5. The process continues until all of the important landmarks are drawn in.
6. Using other colours, the main water supplies are located on the map: i.e. location of water sources/ public taps used by those without household connections to the piped system, sources used during water shortages, spring catchment, tanks and standpipes as well as any wells or rainwater tanks. For piped systems, the pipeline or area of ward covered by this service is drawn in.
7. Again another colour is used to locate sanitation facilities, and problem areas (eg. broken sewers, waste dumps) that might be affecting the sources of water & piped systems.
8. When the map is completed to everyone's satisfaction, it can be redrawn if it becomes too cluttered. This map will be used several times during the course of the training so it is important that it is neat to begin with.
9. Reflection and synthesis of information: Looking at the map, the facilitator asks the group to give their impression on the reliability and quality of the water supply for the various sections of the ward (e.g. going from North-East-South-West), and synthesise their assessment with scores for the points below:

Access: Overall proportion of people with access to piped water system and variation of service in different parts of ward (eg. Access = 1/4-1/2; 1/2-3/4; 3/4-100%; 100%)

Reliability (Water Quantity): Degree to which piped water needs are met during the year by men and women, for different water sources; frequency of service interruption; predictability and influence of irregular delivery.

Water Quality Setting: perceived water quality of different sources identified, protection of water sources, water treatment.

Note: see examples of scoring tables in Annex: Transect Walk. These scores should be tested/adapted at beginning of program.

10. A follow up step is to go to the field to verify the location of the key water and sanitation facilities and problem areas (see section 2.3 "Transect Walk"). It is at this time when the facilitator and sanitarians note (for later use) potential sources of contamination or serious environmental deficiencies that may contribute to poor water quality. The facilitator proposes the organization of this exercise and requests a group of volunteers (men and women) to go along the walk. The group discusses the path the walk will take and decides on a time and meeting place convenient to all parties.
11. Finally, the facilitator asks the group to identify (if any) the location of sites of concern believed to cause contamination of water and or health hazards to the community.

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2.2 Water Use Survey

Purpose: The purpose of this exercise is to determine how the community uses water and what sources are used for each activity.

Proposed Participants: The community at large (men, women and children), the water user's committee (if it exists), relevant NGO and government personnel.

Learning Objective: It is expected that this exercise will reveal water use habits that may be unsafe or contribute to the contamination of the water supply or water used for consumption. Such information is necessary for the design of a health and sanitation education (HSE) program.

Materials:

- ◆ About 20 - 30 medium sized envelopes
- ◆ Photos or drawings of common uses of water: i.e. washing the baby, washing clothes, washing dishes, washing vegetables, cooking, drinking, defecation, water garden, bathing, brushing teeth, etc.
- ◆ Photos or drawings of the water sources that exist in the ward and which are used at different times of year (including during periods of water shortage: i.e. stream, irrigation canal, well, standpipe, household tap, yard tap, etc).
- ◆ beans, seeds or card stock of at least two different colours

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1. At the start of the session a discussion is initiated on what people use water for, and where and how it is stored in the home. Ideally, someone in the community who can draw makes a sketch illustrating each use & storage means.
2. This is followed by a discussion of the different water sources available to the community. Again sketches are made illustrating each. If this is not possible, photos can be taken prior to the meeting.
3. The drawings are then taped to the wall so that the water use habits are across the top and the uses are along the left side in a column. Envelopes are then taped under each picture at the top and to the right making a matrix of rows and columns. See Annex for example. An artist hired by the project can be used to make stylized drawings of water sources and common water use habits.

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4. To start the exercise, all the women are given a handful of coloured beans and the men another. Each participant is asked to go to the chart one at a time to cast a vote by placing a bean in the envelope that shows the water source they themselves regularly use for each activity. It is important to note that if a person is not involved in an activity (eg. if men do not cook) that they should not vote in that column. Only those activities that the person actually does should be considered. If children are present, they can be given a third bean colour and ask to participate in the voting.
5. The voting should be in secret, in other words the matrix is put in a place where no one can see how each person voted, like the hallway. This is so each participant is not influenced by the others.
6. At the end of the exercise, the votes are counted and the results written on the envelopes.
7. A discussion is then initiated on what the group thinks might be an unsafe practice and why. For example: If some people bathe, defecate or dump garbage in the stream and others brush their teeth, wash their dishes, or even drink from the same stream, this is dangerous. See Annex for other possible discussion leads.
8. This leads to a discussion on what the best source is for drinking and why. At this point there is no discussion of boiling water for drinking purposes, or other means of treatment. This comes later in the meeting. Drawings of typical water storage containers are also kept for the discussion of remedial actions.

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2.3 Transect Walk

Purpose:

- 1) To quickly assess the state of water supply and sanitation services in the wards
- 2) To cross-check/ complement some of the information on the community map and water use survey.

Proposed Participants: Interested community members (men, women and children), the water user's committee (if any), relevant NGO and government personnel.

Learning Objective: Make systematic observations while walking through the community from the source(s) of water system(s) to selected delivery points. This exercise will help collect background information necessary for the design of a local WQM program.

Materials:

- ✦ Observation checklists developed for the walk (see examples in Annex)
- ✦ Semi-structured interview guide developed for the walk (see example in Annex)
- ✦ Notepad & pencils
- ✦ 5-10 ready-to-use H₂S test containers.

Timing: 3 x 45 min.

Methods: (should be adapted to the experience/ preference of the facilitator & community context)

1. The exercise is planned at the end of the first meeting (section 2.1 "community Mapping"). A group of men and women representing the ward and interested in water quality monitoring agrees to participate in the walk. At that time, the path of the walk and key places to visit are decided. A suitable time and meeting place are set.
2. The team makes systematic observations while walking from the source(s) of the community water system(s) to selected delivery points.
3. During the walk, the group observes the quality of the water supply and sanitation facilities using the observation checklists. Their observations are discussed with the community members and explanations and insights recorded. Water samples are taken at selected key sites (where contamination is suspected or key public locations). Households (at random) are interviewed along the walk on the availability of water, scope and nature of use, conflicting demands, alternative water sources through out the year and water quality perceptions (see example of semi-structured interview guide in Annex).

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4. After completing the walk, the group returns to a suitable meeting place to discuss and summarize the state of the water systems and sanitation facilities observed, as well as the views of different community members regarding the use and access to services, adequacy and regularity of water supply, and different water sources used through out the year. This information (along with the results of the H₂S testing) will help guide the planning of the local WQM program (sections 3 & 4).
5. The information collected through the first transect walk can be used to refine the observation and interview guides and develop a scoring system of observations and community satisfaction levels which can then be used in designing a more systematic baseline survey for this and other wards.

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3. Monitoring Bacterial Water Quality

3. Monitoring Bacterial Water Quality

3.1 How Water Becomes Contaminated

Lab: Sources of drinking water & how it becomes contaminated

Bacteria in water supplies (What do bacteria look like? How big are they?
How fast they multiply?)

Distribution of bacteria in water (homogeneous vs. heterogeneous)

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Incubation

3.5 Coliplate™ & Colistrip™ Tests

Lab: testing of water samples with Coliplate™ & Colistrips™

Calculating coliform and E.coli concentrations

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H₂S paper strip test

Coliplate™ & Colistrips™ tests

3.7 Cleaning & Disposing of contaminated samples

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3.1 How Water Becomes Contaminated

Lab: Sources of drinking water & how it becomes contaminated
Bacteria in water supplies (What do bacteria look like? How big are they?
How fast they multiply?)
Distribution of bacteria in water (homogeneous vs. heterogeneous)

Purpose: The purpose of this exercise is to show the participants why regular bacterial testing is needed.

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Learning Objective: It is expected that this exercise will show how easily water can become contaminated, that clear water is not necessarily clean water and that disease can be spread by the faecal/ oral route with water as the transport vehicle. Participants will also be introduced to basic characteristics of bacteria and their implications in water quality testing.

Materials:

- ◆ Materials for the H₂S test; i.e. sterile sampling tubes, etc. (See Annex)
- ◆ A wash basin
- ◆ A pitcher of boiled water, covered
- ◆ A flip chart, black board or large paper with masking tape

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1. First, the facilitator initiates a discussion on what water source the group thinks is safe to drink in the city. The facilitator asks the question: "How do you know it's safe?" The answer is: the only way to find out it is safe is to test it.
2. Depending upon the group, a discussion can be initiated on why drinking water may not be safe, linking it to the faecal/oral route of disease transmission and testing for bacteria.

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What Are the Causes of Diarrhoea?

(List all answers on the board or flip chart)

Expected answers: There are unseen organisms that can make you sick if you ingest them. *OR* There are tiny bugs that you can't see that will make you sick if they get into your mouth and you swallow them.

Where Do You Think These Organisms Come From?

(List all answers on the board or flip chart)

Expected answers: From infected faeces. From another person. From human waste.

Microbial versus chemical contamination: Disease-causing bacteria may cause illness within 24 hours or less; chemicals may cause disease such as cancer but it may take many years of consuming this water. This training will deal only with testing for microbial pollution. In other words: testing for the presence of faecal matter in water that indicates that it is unsafe to drink.

3. This leads to a discussion on **how the H₂S test works**. For example, the facilitator explains that there is a simple test to determine if water is contaminated, the H₂S test. Inside the tube is a chemical that will permit bacteria to multiply. If the bacteria that come from faeces are present in the water, the sample tube will turn black within 24-48 hours. If they are not present, the tube will remain clear. Show a tube that is positive and one that is negative. The details of the test will be explained in section 3.3.

4. **Lab Demonstration:** Ask someone from the group to explain what kind of water is definitely safe to drink and why.

Expected answer: Boiled water and bottled water. Boiled water because the organisms were killed and bottled water because is supposed to be tested.

5. Ask someone to put some boiled water in a sample container to the 10 mL mark. Make sure the water is boiled and that no one touches the mouth of the tube.

6. Fill a second sampling tube to the 10 mL mark with bottled drinking water. Make sure the sample container is labelled (with a sample number) and note the type of water sampled on the board with the date and time of sampling. See sample report format in the appendix.

7. Next, the boiled or bottled water is poured over three or more volunteer's hands letting the waste water fall into a basin. Another volunteer takes a 10 mL sample of this water from the basin.

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8. The facilitator says, **“lets go out and find more water to sample.”** One sample tube is given to 5 volunteers. Any place where there is a source of water can be sampled. The location of the sampling point, time of sampling, and sample number are recorded in the participants field notes. The water that can be sampled may include rainwater, ditch water, well water, tap water, water in household storage containers, etc. These tubes should be pre-numbered so that they don't get mixed up.
9. The group returns to the venue and the sample numbers and locations are recorded on the board. The samples are then put in a dark place to incubate overnight.
10. A discussion is then initiated on **basic characteristics of bacteria** and their implications in water testing. Participants are encouraged to volunteer answers to simple questions - allowing the facilitator to gauge the level of knowledge of the group.

What do bacteria look like?

The facilitator invites volunteers to come to the front and draw a bacterium.

Expected outcome:

- bacteria vary in shape (spherical cocci, rod-shaped bacilli, etc...)
- bacteria are very small single celled microorganisms found in our hands, mouth, intestines, clear water, treated water, piped distribution system
- not all bacteria cause disease
- in water testing we are concerned in detecting bacteria that come from faeces of humans and warm blooded animals

How big are they?

Once some drawings are made, the facilitator ask participants to volunteer sizes of the bacteria.

Expected answers:

- bacteria vary in size: 0.5 - 15 Fm (1 Fm = 1/1000 mm)
- single bacteria can be observed only with the aid of a microscope
- in water testing we detect bacteria by allowing them to multiply and grow in the test container: we can see them because we are looking at millions and millions of them in one bottle or tube

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How fast they multiply?

Expected answers:

- bacteria reproduce by dividing: one bacterial cell divides into two cells
- this division can happen at rates up to 100 times faster than with plant cells
- Time for one cell of *Escherichia coli* to divide into two:
 - ◆ 1-2 hours in relatively clean water at 20 deg. C
 - ◆ 20 minutes in friendly environment at 37 deg. C
- One bacterium dividing every 30 min. will produce how many bacteria after 12 hours?
 - ◆ answer: 16, 777, 216 bacteria! (2^{24}).

What is their distribution in a volume of water?

In other words, if there are 1000 bacteria in one litre of water, will every mL that you take out from this volume have one bacterium in it?

Expected answers:

- bacteria are not distributed evenly in water like a cube of sugar dissolved in tea (not homogenous distribution)
- bacteria come in clumps (heterogeneous distribution)
- one single test might not tell the whole story: one volume may contain many bacteria, next volume none, therefore need for regular testing
- coliform densities in distributions system vary - bacteria are not randomly dispersed
- regular testing for presence/ absence of bacteria (frequency) may give better indication of contamination than sporadic tests with quantitative techniques (concentrations or numbers of bacteria per unit volume)

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What indicators & why

What they don't do

Description of two systems: H₂S paper strip & Coliplate™.

Purpose: The purpose of this exercise is to introduce participants to basic principles behind the microbial tests (H₂S paper strip & Coliplate™).

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Learning Objective: At the end of this lesson the participants will understand: how bacterial tests work (need for bacterial growth under optimum conditions, incubation time, incubation temperature); the difference between positive and negative results; difference between qualitative & quantitative testing; and, advantages and limitations of H₂S paper strip & Coliplate™ tests.

Materials:

- Materials for the H₂S test: i.e. sterile sampling tubes, etc. (See Annex)
- Materials for the Coliplate™ and/or Colistrip™ tests
- A wash basin
- A pitcher of boiled water, covered
- A container with contaminated water
- A flip chart, black board or large paper with masking tape.

Timing: 1 x 45 min.

Methods:

1. Referring to the discussions of the preceding section (3.1), the facilitator asks participants to volunteer answers to the following questions:

Given that you cannot see bacteria with the naked eye, how can you test for them?

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Expected answer: you make them multiply in the test container, this means that you need?

- food to grow
- time to grow
- optimum temperature to grow

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How can you make the test tell you if bacteria come from faeces of humans and animals?

Expected answers:

- make a test environment that favours growth of enteric (intestinal) bacteria
 - include a specific indicator to show presence of enteric bacteria (e.g. colour change, fluorescence)
2. The facilitator **reviews rapidly the H₂S test**: Impregnated in the paper strip inside the tube is a chemical that will permit enteric bacteria to multiply quickly. If this happens, H₂S is produced (thus the name of the test) and reacts with another chemical also impregnated in the paper that makes the contents of the tube turn black. More in-depth explanations can be advanced depending on the technical knowledge level of participants.

3 Qualitative & semi-quantitative testing:

- The H₂S test is a qualitative test since it gives YES - NO results (black or clear)
- When the contents turn black, we say that the test results are **positive**
- If there is NO black colour, we say that the test results are **negative**
- If the contents of the test turn black it means that at least one indicator bacteria was present in the test bottle or tube at the beginning of incubation.

- What can be concluded if:

a 100 mL test volume turns **black**? There is at least 1 bacteria in 100 mL.

a 10 mL test volume turns **black**? There is at least 1 bacteria in 10 mL (10 bac/100 mL)

a 1 mL test volume turns **black**? There is at least 1 bacteria in 1 mL (100 bac/ 100 mL)

By varying the sample volume, the H₂S test can be made semi-quantitative

4. The facilitator then introduces the **Coliplate™ test** & demonstrates how to use it (see guide in Annex). The following points are highlighted (these will be reviewed in more depth in section 3.4).

- ◆ Commercial test (ready to use)
- ◆ Quantitative test to measure total coliforms and *E.coli* (common indicator bacteria in WQM)
- ◆ Micro-wells contain special reagents, which make a blue colour appear if coliform bacteria are present and which show fluorescence if the bacteria are *E.coli*.

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- Application is simple, just dispense sample in micro-plate, incubate at 35 deg. C for 24 hours and count micro-wells with blue colour (indicative of coliforms), observe under long wave UV light and count the number of wells that turned blue and fluorescent (indicative of *E. coli*).
- Refer to (MPN) table to determine the Most Probable Number of coliforms and *E. coli* in 100 ml. of sample.

5. Advantages & Limitations:

H₂S test:

- can be produced from first reagents available in-country
- easy to carry out & easy to interpret results
- inexpensive
- prepared paper strips can be stored at room temperature in dark for more than 6 months
- semi-quantitative

Coliplate™ test:

- quantitative test
- specific for *E. coli* (but need UV lamp)
- must be imported & need foreign exchange
- expensive

6. Cost comparisons:

Indonesia (Kromoredjo & Fujioka, 1991):

Colilert test:	US\$6.50	(five-tube tests, no shipping costs included)
P/A + MUG:	US\$1.62	(five-tube test, no cost of preparation included)
H ₂ S test:	US\$0.62	(five-tube test, no cost of preparation included)

Panama (Sánchez & Dutka, 1997):

Coliplate test:	US\$2.70	(US\$4.3 in Jan. 2001) - no shipping costs included
Colistrip test:	US\$1.00	(US\$2.6 in Jan. 2001) - no shipping costs included
H ₂ S test:	US\$ 0.05-0.12	(one tube test, reagents & supplies only, depending on # produced)

7. At the close of the session, two tubes are given to 4-5 volunteers to collect samples at home. These tubes should be pre-numbered so that they don't get mixed up. One of the samples should be boiled water, which is used as a control.

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3.3 The H₂S Test

Lab: basic safety & hygiene in the laboratory

Lab: preparing media & paper strips

100 mL, 20 mL, 10 mL & 1 mL sample tests

Strip preparation versus bottle or tube preparation

Lab: Getting test tubes & bottles ready for use

Negative controls

Purpose: The purpose of this section is to show how a low-cost water quality test works and how it can be prepared.

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Learning objective: At the end of this lesson the participants will understand how the H₂S test works, how to prepare sample bottles and tubes and how to sterilize the containers for use.

Materials:

- Enough sample test tubes with heat resistant lids for about 20-30 samples
- Any type of clean glass bottles with volume of 150 mL to 200 mL with heat resistant lid
- Absorbent paper (non toxic or Whatman filter papers)
- A steamer for sterilization or autoclave (if available)
- Reagents to prepare the H₂S media in liquid form
- Triple beam balance or other scale
- 2-3 pots of at least 500 mL capacity
- Pipettes 5 mL & 1 mL clean & sterile
- Graduated (sterile) cylinders (100 mL)
- A 10 mL pipette graduated in 0.2 mL intervals or a small graduated cylinder
- Marking pencil wax and /or permanent marking pen
- Masking tape
- A ruler graduated in centimeters and millimeters
- Paper towels
- Paper foil (Aluminum)
- Hot air oven to dry papers
- 2 to 3 liters of distilled, bottled or dechlorinated tap water
- 300 g of sugar or table salt (for practicing how to use balance)

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Timing: 3 x 45 min

Methods: (use Module 7: Water Quality Control Techniques in Annex as a guide)

1. The facilitator goes first over **safety and hygiene precautions** while working with microbial testing. (A handout is provided in the Annex for distribution)

Hygiene and Safety Precautions:

The following precautions are essential, not only to *avoid contamination* of the samples and culture medium, but also to avoid the *danger of infection* to staff conducting the microbiological tests,

- ◇ Keep finger nails short; use hair and beard protection if necessary
- ◇ Wash hands with soap and warm water before and after microbiological examinations as well as after use of the toilet.
- ◇ Do not touch your mouth or eyes while working with the test bottles.
- ◇ If you cough or sneeze, cover your mouth. Clean and dry your hands before continuing work with the microbiological preparations or examinations.
- ◇ Keep your work clothes clean.
- ◇ If you feel sick you should not work with the tests.
- ◇ Do not eat, drink or smoke in the working area.

Cleaning the working area:

The working surfaces must be cleaned and disinfected *before and after* each use.

- ◇ Washing the surface with detergent and water removes dust but does not eliminate all germs. It is essential to disinfect the working area. A cheap solution to disinfect all surfaces after they have been washed is to use bleach or chlorine (home disinfectant with 5.25% of available chlorine) mixed with water.

Cleaning solution

Place the bleach (or domestic chlorine) in a clean container and then add the water in the following proportions:

BLEACH		WATER	
5 ml (1 tea spoon)	plus	0.5 liters (2 cups)	
15 ml (1 table spoon)	plus	1.5 liters (6 cups)	

- ◇ The chlorine solution must be left to rest for at least 30 seconds on the surface to be disinfected to ensure this is well done.

In case of accidents or spills: Refer to handout in Annex for instructions

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2. The next step is to **prepare the H₂S media**:

The concentrated medium used in the test is prepared from the following ingredients. These are dissolved by stirring into distilled or non-chlorinated tap water.

H₂S media:

Bacteriological peptone.....	40.0 g
Dipotassium hydrogen phosphate.....	3.0 g
Ferric ammonium citrate	1.50 g
Sodium thiosulphate	2.00 g
Teepol	2.0 mL
Water, distilled or boiled tap.....	100.0 mL

Show participants how to use a balance and allow them to practice weighing reagents (use table salt or sugar for example) and how to use pipettes (use water) before letting them prepare the media.

3. **Prepare the paper strips.** These are the carriers for the H₂S media. The strips can be prepared by either placing them into test tubes or bottles and adding the H₂S media into them - or - by placing the paper strips on sheets of aluminum foil paper and adding the media into them. The last option allows for the preparation of a stock of impregnated papers that can be stored for later use. Both options are outlined below.

- The H₂S test allows to estimate different degrees of contamination by using different volumes of water samples with higher or lower amounts of H₂S media impregnating the paper strips. The following are recommended variations.

Test Sample Volume added to paper strip	Volume of H ₂ S media	Recommended Use
1 ml.	0.5 ml.	untreated waters or waters with suspected contamination
10 mL	0.5 mL	untreated waters
20 ml.	1 ml.	treated waters
100 ml	2.5 ml.	treated waters in well maintained systems

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Preparing paper strips inside test tubes & bottles:

- ◆ Taking non toxic absorbent paper, or Whattman® # 3 filter paper, place a sufficient amount in each container so as to allow the paper to readily absorb the required amount of H_2S culture media. Refer to table above for the volume of media required. Whatever absorbent material is used, it must be cut to a size that will be able to absorb all of the media. This may require some experimentation beforehand.
- ◆ The tubes or bottles can then be loosely capped (do not tighten the caps) and autoclaved for 15 minutes at $115^{\circ}C$. Then, with the caps still lightly loose, dry them in an oven at $55^{\circ}C$.
- ◆ Any type of glass bottles can be used for the 20 and 100 mL tests, as long as they are easy to clean, with a capacity of 50 ml to 200 ml, and with a heat resistant cap.
- ◆ For the 10 mL and 1 mL tests, use 16 x 150 mm test tubes with heat resistant screw caps (or equivalent).
- ◆ Once sterilized, and after they have cooled, the caps should be tightened and the tubes & bottles marked with a line indicating the appropriate volume of sample to be filled (see below). The bottles should be stored in a clean dark place. They can be stored for up to six months.

Preparing stocks of paper strips (dry oven method):

- ◆ Place several pieces of non toxic absorbent paper, or Whattman® # 3 filter paper, on an aluminum foil sheet. The size of the paper should be big enough so as to readily absorb the required amount of H_2S culture media. Refer to table above for the volume of media required. Whatever absorbent material is used, it must be cut to a size that will be able to absorb all of the media.
 - ◆ Place the aluminum sheet with impregnated paper strips in an oven to dry. It is important that the oven is not too hot as the papers may scorch or catch fire. Drying the paper strips at about $70^{\circ}C$ ($160^{\circ}F$) for 20-25 minutes should be sufficient. They should be kept long enough until they are dried but not too long as to burn them brown.
 - ◆ The dried pads or strips are put into a sealed plastic bag, an envelope, or clean glass jar and sealed from humidity for later use. Reagent impregnated paper strips can be stored for several months if kept dry and in the dark.
5. While the paper strips are drying, **the H_2S test procedure** is demonstrated. If the paper strips were prepared inside the test tube or bottle and sterilized, skip the point 6 below and go directly to point 7: "marking of test volumes".

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6. Sterilizing the tubes with paper strips: before going to the field to collect water samples, test tubes and/or bottles with paper strips must be sterilized as follows:

- ◊ Wash the caps and tubes carefully and allow to air dry (use soap, the test tube brush, and rinse thoroughly with clean bottled or boiled water).
- ◊ For each test tube or bottle you are preparing:
 - Cut one dry impregnated paper strip into smaller, thin strips so they can easily fit through the mouth of the test tube or bottle.
 - Carefully place the thin strip pieces inside the tube or bottle.
 - Loosely fit a cap on each tube or bottle
 - Cover each cap with a small piece of tin foil to protect it from the heat (this helps the cap last longer).
 - Place the covered, capped tubes or bottles with the strips in them onto a baking sheet, then into an oven preheated to 110°C.
 - Keep the tubes/ bottles in the oven for 5-6 minutes. Be careful not to burn the paper strips.
 - Carefully remove the tubes/ bottles from the oven and allow them to cool for 10 minutes.
 - Tighten the caps and remove the tin foil.
- ◊ Sterilized tubes/ bottles with paper strips can be stored in a dark place at room temperature for up to 12 months.

7. Marking of test volumes for bottles & test tubes (refer to drawings in Annex handout).

- For tubes which are going to be used to test 10 mL water samples, add 10 mL of water (with the help of a 10 mL pipette or small graduated cylinder) to one empty tube and using a permanent marking pen make a mark on the tube at the maximum height of the added water. Using this mark as a guide, prepare as many tubes as needed with a 10 mL mark line (masking tape can also be used to do the markings).
- For tubes which are going to be used to test 1 mL of water sample, add 9 mL of water to one empty tube and using a permanent marking pen make a mark on the tube at the maximum height of the added water. Then add 1 mL of water to the 9 mL already added and make a mark at the new total level (i.e. the 10 mL level). Using these 9 mL and 10 mL marks as guides prepare as many tubes as needed.
- Similarly, for 20 mL and 100 mL tests, bottles are marked with a line or tape at the appropriate volume mark (use a graduated cylinder to measure the volumes in the guide bottles).

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8. The test tubes and bottles are now ready for use.

Note: Sterilization can also be done by steam (e.g. in a rice steamer) for 30 minutes. If an oven is used it should have a thermometer in it if possible so that the temperature can be monitored. The temperature should not rise too much above 110 deg. C to prevent melting the plastic cap or scorching or burning the paper strip. When cooled they can be stored in a clean place until needed.

9. While the tubes are being sterilized the **concept of a control** is first introduced. It is important that every time samples are collected, a control sample is tested. The control can be either boiled water or bottled water, in other words: water that is known to be suitable for drinking.

10. The water that is used as the control is also used for dilutions i.e. 9mL of dilution water and 1 mL of sample. Every time a 1 mL test is performed, the dilution water should be tested as a control.

NOTE: The facilitator explains that when a control is positive it means that either the tube was accidentally contaminated during sample collection, the sterilization process was not effective or the water used for the control was contaminated. In these situations, the results obtained cannot be considered valid.

11. At the close of the session the facilitator asks for questions. It is important that the concept of the test is well understood as well as the basics of aseptic technique.

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3.4 Sample Collection & Processing

Lab: Water sample collection, tap, well, stream and chlorinated source

Lab: Sample processing

Incubation

Purpose: The purpose of this section is to show how a low-cost water quality test works and how it should be performed.

Proposed Participants: Ward-level participants that will be involved in WQM monitoring, water user's committee, relevant NGO and government personnel.

Learning objective: At the end of this lesson the participants will have learned proper procedures to collect water samples and carry out the H₂S test.

Materials:

- ◇ Enough ready-to-use sample test tubes/bottles with heat resistant lids for about 20-30 samples
- ◇ Marking pencil wax and /or permanent marking pen
- ◇ Masking tape
- ◇ A ruler graduated in centimeters and millimeters
- ◇ Paper towels
- ◇ 1 to 3 liters of boiled or bottled dechlorinated water
- ◇ Incubator
- ◇ Data recording sheets (see Handouts in Annex)

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1. **Choosing the test sample volume:**

- ◇ In situations when bacterial water contamination is not very likely, such as water samples from a properly working system and which includes chlorination, two large sample volumes (20 mL or 100 mL) are used in adequately marked test bottles.
- ◇ For untreated water samples, or waters in which high contamination levels are suspected, two small volumes of sample (10 mL and 1 mL) are used in adequately marked test tubes.

2. **Collecting the samples:**

- The sterilized bottles or test tubes must be opened **only immediately before** collecting the water sample.

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- ◇ If the sample is **non-chlorinated water from tap**, the sample can be collected directly into the test tubes or bottles previously prepared. Turn on the tap and let the water run freely for 30 seconds. Then place the opened H₂S paper strip bottle or test tube under the tap and fill it up to the appropriate volume mark. Close the container, shake it well and incubate as soon as possible.
- ◇ Due to the very small volume to be collected, 1 mL tests require an additional preparation step before the sample can be taken. Use in this case local, boiled tap water. Carefully, through a sterile technique, add the boiled water until it reaches the 9 mL mark BEFORE taking the sample. Then the water sample is carefully added until it reaches the 10 mL mark.
- ◇ If the sample to be tested is **non-chlorinated water from** a well, spring, or storage tank, take the sample using the same container that the home owners normally use to collect the water, rinsing the collecting container at least three times in the water to be collected (**do not put the rinsing water back** into the water source, well or tank). Carefully pour the required volume of test water into the test tubes or bottles utilizing this rinsed utensil. Cap, stir and incubate as soon as possible.
- ◇ To simplify the transfer of sample from the large container to the 10 mL and 1 mL test tubes, pour water from the container into a cup or glass which has been sterilized by having boiled water poured into it and letting it stand for 2 to 3 minutes. Pour out the boiled water, add the sample water and follow the above described procedure.
- ◇ If the water sample comes from a **chlorinated water source**, take the sample in a sterilized bottle (of 120 ml to 150 ml capacity) containing 0.1 ml of 3% sodium thiosulphate (Na₂S₂O₃). This chemical is used to neutralize the chlorine in the water sample. The thiosulphate solution should be added to the sampling bottle prior to sterilizing. After collection of the sample, keep the bottle with water sample in ice or refrigerate it after collecting until such time it can be processed. Iced or refrigerated samples should be tested within six hours. Non-refrigerated samples should be tested as soon as possible. To process these samples, carefully pour the water sample directly into the bottle or test tube that contains the indicator paper strip until it reaches the appropriate mark. Before pouring the water, flame the neck or opening of the sample collection bottle. Cap the paper strip bottle or tube, stir it well and incubate, as soon as possible.

3. **Labelling of tubes and bottles:**

- ◇ The bottles and test tubes must be labeled immediately before each sample is taken to allow them to be properly identified. The label should include the following information:
 - ◇ Identification number.
 - ◇ Origin of the sample (place, type of installation, etc.).
 - ◇ Volume of the sample.
- Date and time the sample was taken.

This information must also be recorded in a field data sheet (see Handout in Annex)

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4. Preservation and incubation of samples.

- ◇ When the samples are collected directly into bottles or test tubes (with paper strips), these samples must be incubated as soon as possible at a temperature between 26° C and 37° C. The preferred temperature is 35° C. A low price incubator, such as the one described in the Annex (Trainee Handouts) can be used. Incubation should continue for a maximum of three days.
- ◇ When the samples are collected in bottles that have been sterilized with sodium thiosulphate (without a paper strip), they must be stored in ice or refrigerated until they can be processed, such as was previously described.

Preparing the Negative ("N") Control:

- ◇ Open a small bottle of non-carbonated water just before pouring it into the test tube.
- ◇ Add 10 mL of bottled non-carbonated water into a tube marked (N) up to the top mark. You may also use water that has been boiled for one minute (then cooled to room temperature) as the control.
- ◇ Immediately cap the tube and place it in the incubator with the other sample tubes. Cover all tubes with a piece of foil paper to protect them from the light of the incubator (see illustration in the Trainee Handouts in Annex).

Notes and Hints:

- ◇ When the sample volumes are only 10 mL or 1 mL, be very careful not to overfill the tubes. If you happen to add a little more your results will still be valid.
- ◇ Timing is important, so plan ahead. The H₂S test should be started right after the sample is taken whenever possible. You should have sterilized and marked the tubes before going to the field to collect your samples.
- ◇ Never use moldy or darkened H₂S paper strips to do your tests.
- ◇ When you take out a new paper strip from its glass jar or plastic bag, make sure you close the containers tightly again as soon as possible. This will help keep the remaining paper strips as dry as possible and allow them to last longer.
- ◇ If the negative control test shows any darkening (a positive result), the tests results of the other samples you processed along with it are invalid. A darkening of the strips in the negative control test tube means that something in the way you are doing the tests is not working. You cannot draw any conclusions from your testing. You should repeat the testing if possible making sure you do not contaminate anything.
- ◇ Make sure that the temperature in the incubator does not rise much above 39° C (higher temperatures may kill the bacteria you are trying to detect.)

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3.5 Coliplate™ & Colistrip™ Tests

Lab: testing of water samples with Coliplate™ & Colistrips™
Calculating coliform and *E.coli* concentrations
Quality Control

Purpose: The purpose of this section is to show how a commercial water quality test for coliform bacteria and *E. coli* works and how it should be performed.

Learning objective: At the end of this lesson the participants will have learned proper procedures to collect water samples and carry out the Coliplate™ & Colistrip™ Tests.

Materials:

- * Coliplate™ & Colistrip™ Tests (6 of each)
- * Marking pencil wax and /or permanent marking pen
- * Masking tape
- * A ruler
- * Paper towels
- * 1 liter of boiled or bottled dechlorinated water
- * A container with contaminated water
- * Incubator
- * MPN tables for interpreting results (see Trainee Handouts in Annex)
- * Data recording sheets (see Trainee Handouts in Annex)

Timing: 1 x 45 min

Methods: (see procedures in Trainees Handouts in Annex)

1. The facilitator first reminds trainees of the use of Total Coliforms and *E.coli* as indicator organisms in WQM.
2. The Coliplate™ & Colistrip™ tests are then demonstrated (see Trainee Handouts in Annex)
3. Trainees are then divided into groups and are given one Coliplate™ & Colistrip™ plates to perform the tests. Plates should be labelled properly before testing a water sample.
4. Samples are incubated & results are measured the following day.
5. Trainees are explained that these tests can be used to confirm H₂S paper strip results and establish *E. coli* densities during the pilot programs (for example on about 10% of the samples taken) in order to convince themselves and water authorities of the validity of H₂S test results. The Trainees are explained the use of the MPN tables for interpreting results. This is reviewed the following day when test plates are actually examined and indicator concentrations are estimated.

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3.6 Recording and Interpreting Results

H₂S paper strip test
Coliplate™ & Colistrips™ tests

Purpose: The purpose of this section is to discuss the results of the sampling done in previous days and to show trainees how this test can be used as an educational tool as well as an early warning sign of a contamination problem.

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Learning objective: At the end of this lesson the participants will have learned how the test is used in a WQM program, and how it can be used to educate the community on why they should disinfect and protect their drinking water.

Materials:

- * The test bottles, tubes & plates of all previous testing
- * Marking pencil wax and /or permanent marking pen
- * Masking tape
- * A ruler graduated in centimeters and millimeters
- * A flip chart, black board or large piece of paper
- * MPN tables for interpreting results (see Trainee Handouts in Annex)
- * Data recording sheets (see Handouts in Annex)

Timing: 2 x 45 min

Methods:

1. After reviewing the topics covered in Day 2, the group examines the results of the tests done on (all) previous days. How the H₂S test reacts on different types of water is noted.

Interpreting and recording the results of the H₂S test:

- Inspect the samples daily for up to three days, or for up to 5 days if the temperature in your incubator is below 30 deg. C (since bacteria take longer to grow at lower temperatures).
- The test is considered positive if it shows any blackening of the indicator paper strip inside the bottle or test tube. In this case, mark a "+" on the record sheet on the day the black colouring first appears. If no black colouring is observed mark a "-". The table at the end of this section can be used as a guide for the interpretation of results. Record the results and your observations in the record sheet (see Trainee Handout in Annex).

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- * Remember: when boiled water is used to prepare the 1 mL test tube, a “*negative control*” must also be incubated with the samples. This is a test tube to which only 10 mL of the boiled water that was used to prepare the samples is added.
- * Note: a **negative result** (no darkening of paper strips) is fine. It means there are no indicator bacteria present in the sampled water. A **positive result** (paper strips turn black) means there are bacteria present.

2. The facilitator then illustrates how the test reacts over time: i.e. after 24 and 48 hours of incubation. For example, a sample that is positive on day one is more contaminated than a sample that becomes positive on day two. In other words, the test can also show the level of contamination by how fast it changes color.

QUESTION: If two 10 mL samples are collected at two different places in the distribution system, and one turns black after 16 hours of incubation and the other after 50 hours of incubation, which one is more contaminated? Which one presents the greater risk to health?

ANSWER: The sample that reacted first, as more organisms cause a reaction in less time. For example: The more people you have to plant the rice, the quicker it gets done.

3. The test can also show the level of contamination by the volume of the sample collected.

QUESTION: If 1 mL and a 10 mL samples are collected and both change color at the same time, which one is more contaminated?

ANSWER: The sample with the smaller volume as there are more organisms per unit volume than in the larger sample.

4. Finally, the test can show degree of contamination by the intensity of the color. For example, if the color of one sample is gray and another is black, which is more contaminated?

ANSWER: The sample that is blacker.

5. If the control is positive what does this mean?

QUESTION: If boiled water causes the H₂S test to turn black, what does this signify?

ANSWER: Either the water was not boiled, was not boiled long enough, the container the boiled water was put in was not clean or the boiled water was mixed with unboiled water.

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Table. H₂S Test - Interpretation of Results

Volume of Sample with a positive result	Amount of bacteria per 100 mL	Observations
1 mL	100 or more indicator bacteria	Probably more than 200 bacteria/100 mL if the blackening takes place very fast and very intensively (less than 18 hrs.).
10 mL	10 or more indicator bacteria	Probably more than 100 bacteria/100 mL if the blackening takes place very fast and very intensively (less than 18 hrs.).
20 mL	5 or more indicator bacteria	Probably more than 50 bacteria/100 mL if the blackening takes place very fast and very intensively (less than 24 hrs.).
100 mL	1 or more indicator bacteria	Probably more than 10 bacteria/100 mL if the blackening takes place very fast and very intensively (less than 24 hrs.).

6. **Interpreting the results of the Coliplate™ & Colistrips™ tests:** The Trainees are explained the use of the MPN tables for interpreting results. They are shown the results of all the plate tests performed the day before - and the fluorescence under UV light of *E.coli* contaminated samples. The Facilitator then records on the board or flip chart the most probable number of indicator bacteria for the different tested samples, asking trainees at random to provide the answers for each plate.

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3.7 Cleaning & Disposing of Contaminated Samples

1. Test containers with bacteriological growth

Any sample that shows any indication of bacteriological growth must be sterilized before final disposal. Sterilize in an autoclave the bottles or test tubes, with their contents and loose lids for 15 minutes at 15 to 20 lbs. of pressure. This process kills any bacteria present in the test containers.

If an autoclave is not available, the containers, and their contents, can be disinfected as follows:

- ❖ Prepare a disinfectant solution in a plastic bucket (such as the one described at the end of the section). Fill the bucket up to half its capacity.
- ❖ Carefully pour into a toilet or a latrine any cloudy sample or any sample whose color has changed, which is evidence of bacteriological growth.
- ❖ Immerse each empty test container in the bucket with the chlorine solution. Leave it there for at least 30 minutes.
- ❖ Afterwards, carefully remove the containers from the bucket and rinse with tap water. Dispose of the empty plastic containers in a recycling bin and, if one is not available, in a garbage pail. If the containers are non disposable bottles or test tubes, wash them with soap and rinse carefully before storing.
- ❖ The following day discard the chlorine solution in the toilet and rinse with a lot of water.
- ❖ Always use protective goggles and rubber gloves when handling any strong chlorine solution.

2. Test or sample collection containers without bacterial growth

If the samples do not show any indication of bacterial growth, they can be disposed in a wash basin, toilet or latrine. Rinse very carefully the wash basin as well as the empty containers. Dispose of the empty plastic containers in a recycling bin and, if one is not available, in a garbage pail. If the containers are non disposable bottles or test tubes, wash them with soap and rinse carefully before storing.

Chlorine disinfectant solution: Add a solution of sodium hypochlorite to a clean container and mix it immediately with water in the following manner: 200 mL of bleach (5% available chlorine) for each litre of water - or - 100 mL of sodium hypochlorite solution (10 to 12% available chlorine) per litre of water.

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4. Setting Up A Local Water Quality Monitoring Program

To protect the health of families in the community, it is important to be able to detect the presence of harmful microorganisms in water before it is consumed. If this can be determined on time, corrective measures can be applied thus preventing diseases and epidemics. It is necessary for this purpose to take regular samples, since water quality can change very rapidly. It is also necessary to identify the most adequate way, in each community, to utilize the information obtained from the water testing to protect water quality and human health.

Activities that a water quality program in each community can implement:

- ❖ Water quality testing at some of the most important supply system locations.
- ❖ Inform health authorities and the community of the test results and the measures that they should take to protect their health.
- ❖ Educate the community about the causes of contamination.
- ❖ Plan and implement protection measures for water supply sources.
- ❖ Install adequate treatment systems.
- ❖ Ensure the proper operation and maintenance of supply and treatment systems.
- ❖ Promote public dialogue and cooperation with the responsible government authorities.

The bacteriological tests learned in this course may be conducted by local technicians in collaboration with water and health authorities, facilitating the availability of information to the community. But information alone is not enough. Households and local organizations are also required to implement measures conducive to the improvement and protection of drinking water quality.

The following sections provide basic elements that should be considered in setting up a pilot community-based WQM program. Also refer to **Module 7: Water Quality Control Techniques in Annex** for further reading.

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4.1 Identifying Sampling Points

Purpose: The purpose of this section is to identify points where water samples should be collected. Sample collection technique and the recording of field data are reviewed.

Proposed Trainees: Community participants that will be involved in WQM monitoring, water user's committee, relevant NGO and government personnel.

Learning objective: By participating in the process of selecting sampling sites it is expected that the participants will begin to understanding the logic behind a WQM program. Also, it is expected that the participants will be able to properly collect water samples and carry out a simple sanitary survey at the completion of this training session.

Materials:

- * The community map showing the community's water supply (see section 2.1).
- * Pens and paper to draw a new map, if necessary.
- * Data sheets (see examples in Annex: Trainee Handouts)
- * Several H₂S ready-made bottles / tubes to take samples.
- * Wooden matches.
- * Report forms or note pads.
- * A container to carry the samples.

Timing: 1 x 45

Methods:

1. This is the first step in establishing a local WQM program. Before reaching this point, selected community members were trained on the preparation and use of the H₂S test (section 3), the community decided to participate in a community WQM program, and the mapping and transect walk exercises in the community were carried out (see section 2).
2. The facilitator can begin by going over the data sheets. The form(s) used to record the sampling points and results of testing is put on the wall and explained.
3. Then the tubes of samples collected in the transect walk (section 2.3) are examined and the results recorded with a plus for positive and a minus for negative.
4. The results of the earlier sampling done at home by volunteers is also recorded by those who collected the samples. In this case, the person who collected the sample explains the conditions where the sample was collected and whether or not they expect the results to be positive or negative.

(Note: Depending when these samples were taken the tubes may still be negative. Black colouring often does not start to take place until after about 18 hours of incubation.)

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5. The facilitator then initiates a discussion on the results of each of the tests.

QUESTION: Which water is safe to drink and why? Which water is not safe to drink and why?

6. This leads to a discussion that only after the water has been tested can we tell that it is safe. Therefore samples must be collected from all those places where the water is most likely to become contaminated.

7. The community map (section 2.1) is used to identify the sampling points for the locality.

8. The facilitator initiates a discussion on where the group thinks samples should be collected and why. Some points the facilitator should keep in mind during the discussion:

- * Identifying the sites where samples should be taken, and how frequently, will depend on: the available resources, the number of persons that may become ill through contamination of a specific site, the type of water source and water supply system, and the quality of the water.
- * Properly operating treatment plants, producing good quality water, do not need to be sampled as often as those with operational problems.
- * Water pipes (water distribution networks) that regularly lose pressure (periods without water) must be sampled at more sites and more often than those which always maintain pressure. When "water supply is cut" the pipe can become contaminated with waste water from outside the system.
- * Certain points in the pipes/aqueducts, or the community are more important than others and require more frequent samples being taken. These are places where there are:
 - many children (such as schools, kindergartens, etc);
 - old people (such as retirement homes);
 - sick people (such as clinics and hospitals).
- * Samples should be taken in public wells at least twice a year, and one of these samples must be taken during the rainy season. If the well supplies many families, the frequency of the sampling must be higher, particularly if the individuals are not used to disinfecting their water.
- * In homes with their own wells or where drinking water is kept in buckets, containers or water tanks, samples can be taken just once a year, on a random basis, if the community considers it appropriate and if sufficient resources are available.

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- When water is sold to individual homes from water tank trucks, a sample per truck must be taken at least once a week or once a month, depending on available resources (One should not pay for contaminated water!).
- * It should be kept in mind that, many times, a water pipe system is not sufficient to supply the whole community, either because there is not enough water during the dry season, or because the pipes never reached certain districts of the community. An effort must be made to know the community and the use of the different water sources available (see section 2.2). In general, those not served by a water distribution system (pipes) are also those that consume the most contaminated water and the ones that have the least amount of resources to combat diseases.
9. The sampling points are marked on the map with numbers to identify them. Begin by marking critical public sites (e.g. public wells, stand pipes, health clinics, water trucks, school yard water fountains). These numbers for public sampling points should be kept fixed. A fixed numbering system prevents errors in the event samples are mixed up: i.e. each sampling point is identified by the same number each time. This also allows the water authority and the sanitarians responsible for the community to identify changes in water quality over time, one of the objectives of the WQM program. If private homes are sampled (e.g. as part of hygiene awareness raising), the numbering system may reflect this (e.g. H1, H2, stand for Household 1, Household 2,)
 10. After it is agreed that those points selected should be the sampling points for the ward, the participants accompany the ward WQM technicians to collect samples.
 11. At each sampling station, the sampling technique is demonstrated by the ward trainees. Taps can be flamed, but it is not necessary, as it is the quality of the water that is being brought to the home that is important. A sample can even be taken directly from the container of a ward member collecting water at the sampling site (but this is recorded in the data sheet).

Notes:

- All public wells, the spring catchment and all distribution points including water storage tanks should be sampled. House connections need not be sampled unless done to show that a poor plumbing connection or an unclean household storage container can contribute to a contamination problem.
 - If it is the quality of the water in the piped system that is being checked, then the taps should be flamed. If flaming is done, it is not necessary to use a Bunsen burner. Three or four wooden matches held together are ignited and held under the tap until the water droplets at the mouth of the tap turn to vapor.
12. It is important to note that the sample collector's hands should not touch the mouth of the sample container, nor should he/she put the cap on an unclean surface while taking the sample. A trained sanitarian will be able to demonstrate the correct procedure for opening and holding the container and cap.
 13. Before or after collecting the sample (preferably before), information on the location and number of the sampling station, date and time of day, etc. are recorded.

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14. An inspection (called a sanitary survey) is then done to see if the tank is clean, if there are sources of pollution nearby, or if there are any components in need of repair or replacement - like taps, tank covers, etc. The facilitator and the trainees keep their own notes. This permits the facilitator to see who recognizes what. The various observations are discussed later on in the session and if necessary the sites are revisited.
15. Also, if the water is noticeably turbid this is recorded as well. Checking turbidity can be done by filling a white basin with water. If the water is turbid it will appear milky, brown or colored in the basin.
16. The sample number is marked on the sample container and the group moves on to the next sampling site. At each site a different ward trainee does the sampling with the facilitator looking on to see if there is any accidental contamination; for example, the cap of the sampling tube placed in an unclean place, etc.
17. Upon returning to the training venue, the facilitator initiates a discussion on which sampling points the group think are the most contaminated. What problems were observed like sources of pollution and maintenance problems. What does the group think can be done about these problems?
18. The tubes of the day's sampling are then left to incubate overnight.

Notes:

- * Looking for potential sources of pollution, construction defects or other problems that can effect water quality is called a sanitary survey.
- * It may be necessary to equip the trainees with flashlights so that they can see in dark places properly.
- * Turbidity is an indicator that either surface water and/or pollutants are entering the system.

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4.2 Roles and Responsibilities

Purpose: The purpose of this section is to define a local WQM program, develop an action plan and determine who will implement it.

Proposed Trainees: Community participants that will be involved in WQM monitoring, water user's committee, relevant NGO and government personnel. NOTE: This exercise cannot be done until technical training on the H₂S test (section 3) has been completed.

Learning objective: To gain skills in planning and program development.

Materials: A flip chart, black board or large piece of paper with masking tape.

Timing: 2 x 45 min.

Methods:

The questions of WHO, WHAT, WHERE, WHY AND WHEN are addressed in this training session. This assumes that the question of HOW has been dealt with in previous sessions.

1. The first step is to establish the reason for a local WQM program.

The H₂S test allows you to monitor the quality of the water that people drink through out the year (times of water shortages included) and can prevent illness by alerting you that a contamination problem exists. If the test shows the level of contamination has increased, then this indicates that maintenance or health preventive action is needed. As long as the system is operational, changes are likely to take place and it is the job of the WQM program to detect those changes and take action to avoid serious problems. The program can also be an ideal venue to promote health and hygiene measures in the community.

2. The facilitator then initiates a discussion on what tasks need to be done to monitor the community's water supply system, and writes this down on a flip chart or board.

For example, people are needed for the following tasks:

- sterilize sample containers
- collect samples
- carry out a sanitary survey
- keep track of accounts

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- * keep records and send in reports
- * inform responsible authorities on the monitoring results
- * lead community education sessions and carry out demonstrations
- * recommend and mobilize the community for remedial action or maintenance work
- * supervise the work
- * other duties (list)

NOTE: Many of the above mentioned tasks can be done by the same person. It is recommended that at least a couple of ward members be trained on these tasks to provide support and/ or serve as a back up.

3. Once the jobs have been defined, who does each task is the next question. The people chosen to implement the program should only be those who have the time and the interest to do a good job. The important thing to keep in mind is that the person(s) chosen to run the program is responsible for the environmental health of the community. Just as you wouldn't want to be looked after by an incompetent Doctor, you don't want to entrust the job of monitoring water quality to someone who is really not interested.
4. The next topic of discussion is to define who all the institutional actors are, what their existing roles are and what role they intend to play in the new program.
5. The discussion can begin by listing who should be involved (from being informed about the program to being actively involved in its implementation).

Drinking water quality depends on a great number of factors related to the environment, the type of construction, the operation and maintenance of the supply and treatment system, current sanitary measures and hygiene practices. In turn, the quality of the drinking water has a direct impact over the health of the individual and the community. For these reasons, there are several persons in each community that, one way or another, share responsibilities and interest in ensuring a good water quality and routine monitoring. This may include:

- * organizations or committees related to the water supply and sanitary system (such as a Water Administration Board or Committee);
- * individuals or community groups working on health issues (health committee, health workers, nurses, doctors, and other similar groups);
- * representatives of the local, municipal or state authorities responsible for the water supply system and for water quality;
- * respected individuals (leaders) in the community, such as school teachers, religious leaders, etc.

It is necessary to discuss with all of them, as a group or individually, the objectives, benefits and scope of the program. Their support and advice should be requested when planning and implementing activities designed to improve and protect drinking water quality and community health.

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6. Then, the role of key institutional actors is discussed and key points advanced by the participants are written on the board or flip chart. For example:

a) Local ward WQM committee (& collaborating NGOs during pilot program)

Suggestion:

- * To provide technical support and advise to the ward WQM technician(s)
- * To assist the technician(s) in fulfilling reporting requirements
- * To investigate reported water quality problems and seek corrective action
- * To keep financial records and provide the WQM technician(s) with the materials they need to carry out the program.
- * To collect, collate and analyze WQ data generated by the program within their jurisdiction and send this to the responsible authority
- * To inform the community about the progress/ achievements of the local WQM program

b) Municipal water authority

Suggestion:

- * To provide overall support, program direction, and formulate the operating budget of pilot programs
- * To make the H₂S media and provide this and sample tubes to each participating ward
- * To provide follow-up training to ward technician(s) about the technical aspects of WQ control.
- * To assist in selecting the wards for the "pilot program"
- * To act as facilitators in initiating WQM programs in these wards
- * To provide technical and administrative support to the ward WQM committee
- * To monitor each WQM program and document progress in upgrading their systems
- * To collect, collate and analyze data from each ward and send this information to the responsible authority
- * To provide technical support in the planning and implementation of corrective actions
- * To analyse and document the overall progress of the pilot programs and information generated to assess cost effectiveness and technical viability of the programs

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c) Central planning, water & /or health authority (based on the outcomes of the pilot programs)

Suggestion:

- * To provide overall support, program direction, and formulate the operating budget of an expanded program
- * To define the type of information to be collected from the wards
- * To collate and analyze statistics collected and link this information to health indicators
- * To report on progress in the provision of potable water to the city
- * To report on the status of potable water supplies in the city
- * To analyze the data collected to determine what programs and projects are the most cost effective and technically sound in achieving program goals
- * To petition for funds to expand the program to those communities not yet served
- * To develop new program directions based on the information at hand.

7. To focus the discussion, the facilitator may pose the following question to the group and record the answers in a sheet of paper in the wall:

QUESTION: Why is it necessary to Monitor Water Quality?

There are several answers to this question:

- * To determine changes in water quality that may indicate contamination problems
- * To protect the health of the community by advising people to boil or treat their water when tests indicate a contamination problem.
- * To educate the community on the importance of protecting their water resources from contamination & the need for proper hygiene and sanitation practices
- * To determine where a contamination problem may be located
- * To determine how well sanitary protection measures are working
- * To initiate maintenance when indicated.

8. The next step is to define the frequency that samples are to be collected. It is suggested that initially sampling be done once or twice a month (depending on key sampling points selected). If after 6 months the results are consistent and the system is being well maintained, then the sampling can be done less frequently.

9. The results of this session are used for the next section "Planning the Program."

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4.3 Planning the Program

Purpose: The purpose of this section is to define an action plan and establish a local WQM program.

Proposed Trainees: Community participants that will be involved in WQM monitoring, water user's committee, relevant NGO and government personnel. NOTE: This exercise cannot be done until technical training on the H₂S test (section 3) has been completed.

Learning objective: To gain skills in program planning and implementation and establish a local WQM program.

Materials: A flip chart, black board or large piece of paper with masking tape.

Timing: 2 x 45 min.

Methods:

1. The facilitator asks those present to define their present roles concerning the pilot WQM program and the operation and maintenance of the ward water supply and sanitation systems. Each job is listed and the duties summarized under the job title. If possible, an organogram is drawn showing who reports and/or should report to who.
2. Next, the facilitator asks who can and should carry out the duties defined in section 4.2 "Roles and Responsibilities". The names of these people are added to the organogram.
3. The facilitator then focuses on finances and resources. In these discussions the following issues must be examined:
 - * needs (location, equipment, inputs, transport)
 - * how much will it cost to carry out such a program (look at transport, the cost per sample, fees, etc...)
 - * available resources
 - * limitations (what can and cannot be done)

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4. Finally, the facilitator summarizes the program: i.e. what the community will do, what human and material resources are needed and how much it will cost.
5. This and other relevant program information is written down and becomes the ward's implementation plan. Once everyone agrees on what will be done and by whom, the resulting document is signed by all concerned parties and the person(s) who will be responsible for routine sampling and testing, inspecting the system and reporting the findings. In order to sustain this input, the person(s) selected should be given an appropriate incentive allowance, the amount to be determined by the program organizers.
6. The agreement is kept on file as a "memorandum of understanding." This information will be used to later evaluate the program: i.e. how well the community and government agencies have complied with their role as defined by the agreement.

Note: The "memorandum of understanding" describes the program so that there is no confusion about expected activities, roles and responsibilities.

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4.4 Reporting Requirements

Purpose: The purpose of this section is to establish a reporting system..

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Learning objective: To learn how to develop and implement a records system and reporting procedures

Materials: Reporting forms; flip chart, black board or large paper and masking tape.

Timing: 2 x 45 min.

Method:

1. As part of this exercise, all records and reports (if any) that are already part of the ward's environmental monitoring program should be reviewed and discussed.
2. It is important that the trainees understand that records must be kept.

QUESTION: In the context of WQM, ask the participants what type of data should be collected and why. Record the answers on the flip chart.

Expected answers:

- 1) the results of all testing, to show the status of water quality in the community;
 - 2) the results of sanitary inspections, to show where the problems are located;
 - 3) a description of repairs or record of maintenance undertaken, to determine what breaks, how often and to plan for the future; and,
 - 4) the cost of such remedial action, to be able to budget and account for the money spent.
3. The format for recording the data is discussed as well as the frequency of reporting and to whom this information is forwarded to. It is most important to inform the community why the information is needed and how it will be used. Collecting data for the sake of collecting it can be a disincentive. All parties have to do their job and the information must flow both up and down the chain, if the program is to be effective.

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4. It is also important that whatever reporting format is used, that the data collected, including descriptions of any action taken, can be linked to the water supply in the ward where the samples were taken.

Note: It makes no sense to report that (x) number of samples were collected and (y) number were found positive and (z) number were found negative. The information has to be presented in a way that correlations can be made. In other words, which type of water supply systems have the best quality water and where are they located? Is there a correlation between an active O&M program and good water quality? Has the WQM program made a difference in improving O&M and water quality?

With the appropriate data, answers to questions like these can be obtained. They are not just guesses or anecdotal information, but are backed up by relevant statistics.

5. Examples of sample collection and H₂S test report forms are found in the Annex.
6. Finally, the facilitator should take the opportunity to discuss the reporting procedures for cases where contamination is found. He/ she should begin by stressing the importance of maintaining confidentiality of the water test results and the observations that take place during visits to different private homes and public facilities, except for the purpose of informing the responsible persons or institutions. These may be the family head (in the case of a home visit), the staff operating the treatment plant and the corresponding water authority (if the samples have been collected at the treatment plant or the piped system), the director or the teacher of the school (if the samples were taken in the school), and community health authorities when the samples have been taken in the piped system or in a public place.

Discretion and respect are two required elements to do a good job

It is important to keep in mind that if discretion and confidentiality are not kept, the trust of the community may be lost. If good relations are maintained with the persons visited at their homes or public facilities, this will allow to discuss with them any problem that may arise and the preventive or corrective measures required to improve the water quality.

7. This leads to a discussion of **Notification of Results and Corrective Measures**. The following points may be used as a guide.

What to do when a sample shows positive results (i.e. contamination)?

If the sample taken in a family home indicates the water to be contaminated, the following steps should be taken:

- Visit the home personally or inform the owner by phone and, if the sample comes from a water tap connected to the water distribution network, inform as well the person responsible for the operation and maintenance of that network.

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- ❖ Take another sample from a neighboring location using the same water source (for example, from the tap of a neighboring house connected to the same water supply system or the same well, or the same water tank or public tap).
- ❖ If the second sample also shows positive results, inform all neighboring homes, and the person responsible for the water source in question, of the degree of contamination found.
- ❖ During the follow up visit, examine the possible water contamination sources during storage or handling, discuss the measures that may be taken and recommend that, until the situation has been corrected, water should be boiled or disinfected before drinking.
- ❖ Document the results of the bacteriological tests and their follow up.

If a **sample from a public facility or system** indicates the water to be contaminated, the following steps should be taken:

- ❖ Inform immediately the persons responsible for the facility and the health representative of the community (or the person with whom the work has previously been coordinated).
- ❖ Visit again, accompanied by the person responsible for the system, the location or locations affected and examine potential contamination sources, providing the corrective measures that the system operators or responsible persons can adopt (for example, in a chlorination plant, verify if the chlorination equipment is working properly and make the necessary adjustments: in the case of a storage tank, clean and disinfect it).
- ❖ Collect additional samples and repeat the bacteriological test. If possible send a sample to an accredited laboratory to determine, with their help, the level of contamination.
- ❖ Ask the person responsible for community health to recommend that all users of the supply system (or water source) boil or disinfect the water before drinking, until the situation has been corrected. This step is a must if the second set of samples confirm the contamination of the system.
- ❖ When informed that the problem in the system has been corrected, collect a new sample to verify that there is no more contamination.
- ❖ Document the results of the bacteriological tests and their follow up.

The above steps are given only as a guide. They will have to be modified and adopted in collaboration with water and health authorities, as well as local authorities, to meet local community requirements. What is most important is to have these steps well defined before starting the program.

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5. Some Remedial Actions

Use of regional support and resources

Available resources will determine the best way a problem can be handled. Some solutions may require outside resources (materials, funds or technical support). For example, if the chlorinating mechanism of a treatment plant breaks down, help may be necessary to have it repaired.

There are however several activities that the community may undertake, jointly or in support of an outside intervention. It is necessary for this purpose to use the available resources of the region. It should be kept in mind that, under any circumstances, the most important objective is to prevent water contamination in particular and environmental contamination in general. For example, a group of family members, of neighbours, or the whole community, may get organized to implement activities such as:

Use water quality monitoring as an opportunity to work in contamination prevention

- ❖ Protect and keep clean the area around the water source.
- ❖ Construct a well built latrine and keep it clean and ventilated.

If the people in your community use water from a well that is not well protected, or from the river, or when the water distribution system is contaminated, promote family health protection practices such as:

- ❖ Filter the water through home filters made of fine sand and small gravel (and boil or disinfect the filtered water).
- ❖ Add chlorine to the water in the well or storage tank (consult with the health representative of your community).
- ❖ Boil drinking water.
- ❖ Always keep the storage containers for drinking water covered and protected, and out of reach of animals.

All these subjects are covered in more detail in the respective Modules of this series

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5.1 Filtration

Slow sand filters were one of the first techniques used for treating water. Filters allowed the removal of most impurities suspended in the water and helped prevent the spread of diseases such as cholera.

During filtration, impurities are retained in the sand as the water drains through. They become trapped or adhere to sand particles. However, other purification mechanisms are also involved. With the use of the filter, a microbiological layer forms at the top surface. This layer is composed of many different types of microorganisms that feed on the organic matter carried by the water being treated. The layer of microorganisms that formed contributes to clean the water from impurities as these become absorbed into the layer. Given its beneficial role in treatment, it is necessary to keep at all times sufficient water covering the filter surface so as to maintain the biological layer immersed and ensure the survival and reproduction of the microorganisms it contains.


If the water to be filtered contains pathogenic (disease-causing) organisms, it must be disinfected with chlorine before drinking. In this case, the water must be disinfected **AFTER** filtration otherwise the chlorine would also damage the biological layer of the filter.

In the Annex we describe the building procedures for an inexpensive filter design for home use recommended by Davnor Water Treatment Technologies Ltd. of Canada. A diagram of the filter is shown in the next page.

A recent study conducted by the National Water Research Institute of Canada (NWRI), with support from IDRC, demonstrated that these home filters are capable of removing from the water not just bacteria but also parasites and toxic substances. This study was conducted utilizing waters with toxic concentrations 10 to 100 times larger than those encountered through normal contamination. The removal rates obtained were 83% for total heterotrophic bacteria populations; 100% for *Giardia* cysts; 99.98% for *Cryptosporidium* cysts; and of 50 to 90% for organic and inorganic toxic substances⁶.

⁶ Palmateer et al., (1999). Toxicant and Parasitic Challenge of Manz Intermittent Slow Sand Filter. *Environ Toxicol* 14:217-225.

Manz Intermittent Slow Sand Filter

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5.2 Disinfection

Disinfection of water refers to the destruction of germs to safe levels, making the water potable (drinkable) with no risk of transmitting diarrheal disease when consumed. Disinfection should not be confused with sterilization which refers to the total destruction of microorganisms in the water. Sterilization is much harder to achieve, requires more energy and is more expensive. For drinking purposes, water need not be sterilized, but it should be disinfected. The two most common forms of disinfection are:

- **Boiling:** Boiling is the oldest, safest and most effective method of disinfecting water. It works on the principle that microorganisms cannot tolerate the high temperatures and bacterial cells rupture and proteins are denatured. However, boiling water requires important amounts of energy. It takes one kilogram of firewood to bring one liter of water to boil for a minute. The amount of time that has been traditionally recommended for boiling water varies from 5 to 10 minutes. Such a long period is not only unnecessary but makes this option economically and environmentally unsustainable. Full boil for one minute is sufficient to kill any disease causing bacteria in water and inactivate Giardia cysts.
- **Chlorination:** This is the method used in most cities around the world to disinfect water. Although the taste of chlorine may be disagreeable to those who are not used to it, it is harmless and leaves a disinfectant residual in the water for added protection.

Disinfection Exercise

Purpose:

The purpose of this module is to demonstrate how to clean and disinfect a storage container, cistern, tank truck, well and piped water supply system.

Proposed Trainees:

Ward-level participants that will be involved in WQM monitoring, water user's committee, relevant NGO and government personnel.

Learning objective:

The participants will be able to understand the need for disinfection and the procedures for cleaning and disinfecting various components of a water supply system.

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Materials:

- A supply of calcium hypochlorite (bleaching powder dissolved in clean water-see handling recommendations & preparation below), and sodium hypochlorite (laundry bleach).
- A stiff brush for cleaning the inner walls of tanks.
- A plastic bucket.
- 10-12 H₂S test tubes/bottles (ready to use).
- An empty (but clean) 600ml clear plastic bottle (with cover).
- A 5-10 liter sample of known contaminated water collected from the community's water supply (water should however be as clear as possible - i.e. of low turbidity).
- A pan used to boil water and a source of heat.

Timing: 3 x 45 min.

Method:

1. The first thing to cover is the concept of disinfection. In section 3, drinking water and boiled water were tested with the H₂S test and the results compared with the tests done with contaminated water. This section should have established the fact that there are undesirable organisms in the water and that these are killed when the water is treated. If the water is not treated before drinking, the consumer's health is at risk. In other words, contaminated water can be made safe to drink if it is boiled. This is called disinfection.
2. The facilitator mentions that there are other methods for disinfecting water (making it safe for drinking purposes). Three common methods are covered in this exercise:
 - 1) by heat. (boiling)
 - 2) by chemical means (chlorination); and.
 - 3) by solar energy.
3. To show the three methods of disinfection, an experiment is done.
 - H₂S testing of contaminated water: First, 10 ml and 1 ml samples of contaminated water are taken directly from the container it was collected in. They are labeled and put aside.
 - For solar disinfection: The clear plastic bottles is filled with contaminated water and a 10 ml sample of water in the bottle is taken. The bottle is then put in the sun for a period of at least 4 hours. After 4 hours of exposure another 10 ml sample is taken for H₂S testing. It is important that the bottle receives direct sunlight for the entire period. The bottle should be placed on a white or reflective surface so that the solar energy can bounce back into the container. The time of day recommended for exposure is between the hours of 10:00 AM and 2:00 PM.

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- For boiling: Half of the contaminated water in the original collection container is removed and boiled for ONE minute. When cooled, another 10ml sample is taken for H₂S testing.
- For chemical disinfection: Two drinking glasses are filled with contaminated water. In each glass are added twenty (20) drops of the bleaching powder solution (see how to prepare solution below) OR five (5) drops of laundry bleach. After about 20-30 minutes of contact time, 10ml samples of the chlorinated water are taken from each glass for H₂S testing. They are labeled and incubated.

Important Notes:

- Commercial laundry bleach solutions (primary sodium hypochlorite) and commercial bleaching powders (primarily calcium hypochlorite) used for laundering often contain impurities or additives to improve cleaning. They are used here for purposes of demonstration only. These disinfectants can be used to clean storage containers and tanks (with good subsequent rinsing), but pure reagents for drinking water disinfection should be used to treat water.
- A bleaching powder solution for this exercise can be prepared as follows:
 - Use fresh (less than 6 months old) bleaching powder (65% available chlorine).
 - Add one liter of clean water into a clean container.
 - Add into the water four small spoonfuls (16 grams) of bleaching powder.
 - Mix well, cover container and let stand for 1 hour before use. This solution should contain about 1% available chlorine.
- It must be explained that bleaching powder is a strong oxidizing agent and must be treated with care: i.e. it is dangerous to handle it directly. It can also damage clothing.
- Bleaching powder should always be added to the water, not the water added to the powder.
- Sodium hypochlorite, or bleach used for laundry purposes can also be used for the exercise. This solution usually contains about 5% of available chlorine.
- The amounts of chlorine solution being added to the glasses with contaminated water are much higher than would normally be added to drinking water. This is done only for demonstration purposes only, given that we do not know how contaminated the water is. *This can be a good opportunity for the facilitator to remark that if the dirty water had been filtered through the slow sand filter described in section 5.1, we would need less than one (1) drop in each glass to make the water safe to drink.*

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- If possible, bottled water with a chlorine odor and taste is also purchased for demonstration purposes and sampled for H₂S testing.
- The samples are allowed to incubate overnight and are examined the next day.

RESULTS:

The results should show that those samples treated by boiling, solar radiation, and chlorine have not changed color (did not turn back), indicating that the bacteria that cause the test to react positive have been killed.

The facilitator asks one or several of the participants to explain the results of the above experiment to the group. The discussion continues until everyone is clear about what is meant by disinfection.

QUESTION: How do you know when your water supply system (well or piped system) needs disinfection? Which method do you think is used and why?

ANSWER:

- When the H₂S test shows that the water in the system is more contaminated than at the source.
- When foreign matter or algae is observed in the well or storage tanks.
- After any maintenance or construction activity.
- When the H₂S test reacts positive within 24 hours with a 1 ml sample.

Steps involved in disinfecting various components of a water supply system.

The following are possible procedures for disinfecting various components of a water supply system. They are presented here as examples only. Before recommending them to users, health and water authorities should be consulted and standard procedures defined. This consultation and coordination is especially necessary for the disinfection of public wells, water trucks and sections of the piped distribution system. Due to the nature of the chemicals used to disinfect a water supply, supervision by an experienced sanitarian is necessary when disinfecting public facilities.

Shock Chlorination: This refers to the disinfection of components of the system (not of the drinking water itself). Disinfection is achieved by adding a specified amount of chlorine solution, leaving the chlorine solution in contact with the component being disinfected, draining the chlorine solution, and rinsing the component before using it again. This is called "Shocking" the water component (e.g. *shocking the well, shocking the cistern, shocking the barrel, or shocking the water main*).

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Water Barrels: The size of the water barrel will determine the amount of chlorine required to disinfect the barrel. Suggested amounts are:

Volume of barrel	Amount of bleach (5% available chlorine) added
1 L	1 drop
1 Gallon	2 drops
20 L	1 teaspoon
1000 L	2 tablespoons

Steps:

- empty barrel and clean
- fill with clean water
- add required amount of chlorine solution
- let stand for 30 minutes
- empty barrel and rinse thoroughly
- refill barrel with clean water, disinfect and close lead.

Water Cisterns: The size of the water cistern will determine the amount of chlorine required to disinfect the barrel. Suggested amounts are:

Volume of cistern	Amount of beach (5% available chlorine) added
4,000 L	4 L
1,000 Gallons	1 Gallon
8,000 L	8 L
2,000 Gallons	2 Gallons

Steps:

- Before adding chlorine to cistern, its size must be determined.
- Disinfection procedures should be done for all newly installed cisterns or when bacterial tests show high levels of contamination.
- All cisterns should be disinfected at least once every year.
- Before beginning any disinfection procedures, cut and remove any vegetation or debris from around the cistern. Inspect the cistern lid and fill cover for cracks.

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- Before starting, make provisions for storing enough water to meet household needs for 12 HOURS.
- Shut - off any by-pass the hot water tank and any filters that may be connected to the water system.
- Drain cistern. If possible (and with supervision/assistance) clean the interior walls of the cistern by scrubbing or using a hose to remove foreign material from the walls. Using a shovel or pail, remove all debris from the cistern. Inspect the interior walls for cracks and water leaking into or out of the cistern.
- Refill the cistern with fresh chlorinated water.
- Pour chlorine solution into the cistern as calculated from the table.
- Turn on the cold water faucet at the kitchen sink. Let the cold water run until a chlorine odor is present. Shut off the kitchen faucet and repeat this procedure for every faucet in the house or yard.
- Leave the chlorinated water in the household water system for at least 8 HOURS.
- Connect a garden hose to the outside tap of the house. Place the other end in the drain or area of little or no vegetation. Turn the outside tap and all faucets and let the water run until the chlorine odor is gone.
- Re-fill cistern with fresh chlorinated water.
- After on week, re-sample the tap water for bacterial content (do H₂S test).

Water Trucks: The size of the water tank truck will determine the amount of chlorine required to disinfect the barrel. Suggested amounts are:

Volume of tank truck	Amount of bleach (5% available chlorine) added
2,000 L & Smaller	2 L
500 Gallons & smaller	1/2 Gallon
4,000 L	4 L
1,000 Gallons	1 Gallons
6,000 L	6 L
1,500 Gallons	1.5 Gallons

Steps:

- Disinfection procedures should be done for all newly purchased water trucks or when bacterial contamination is detected in the water.
- All water trucks should be disinfection procedures, inspect all delivery hoses for tares and nozzles for rust.
- Drain water tank.
- Using a flashlight, inspect as much of the tank interior as possible. If rust or debris can be seen, report it to corresponding authority. Any debris must be removed and rust sanded off and re-coated.

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- Re-fill tank with fresh chlorinated water.
- Add appropriate amount of chlorine solution to tank.
- Flush all delivery hoses. Shut off after chlorine odor is present.
- Have the tank truck operator drive the vehicle around for a few minutes (to mix contents). Let stand for THIRTY (30) MINUTES.
- Drain the water tank through the delivery hoses.
- Rinse the water tank with fresh chlorinated water and flush again the delivery hoses.
- Re-fill with fresh chlorinated water.
- Re-sample and perform bacterial testing.

Piped system: This activity should be carried out with the supervision of the water authority responsible for the water supply system.

Steps:

- Before starting, inform all household attached to the portion of system to be flushed to make provisions for storing enough water to meet household needs for 12 HOURS.
- Then all of the holding tanks are scrubbed with a stiff brush to get rid of algae and foreign matter. Use a mild chlorine solution. Do not use soap.
- The tanks are then flushed out to get rid of the debris and organic matter.
- The chlorine compound solution is then introduced at the source.
- Once the smell of chlorine is noted in the water coming out of the taps, the valves are closed and the chlorinated water is allowed to stay in the system overnight.
- The next morning the entire system is flushed.
- After the smell of chlorine has dissipated, samples are taken. The results should show that the water in the system is of the same quality as the water at the source.
- If not, a further investigation is needed to identify the source of contamination, and appropriate corrective action is taken.

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Wells: Cleaning and disinfection of wells is more difficult. Regular monitoring, household disinfection and proper storage of water are essential steps to guarantee their portability.

Steps:

- The users of the well should be told to make arrangements for another source of water while the well is being disinfected.
- First the concrete lined walls of the well must be scrubbed using a weak chlorine solution (e.g. 0.5% available chlorine).
- Next, the water must be removed from the well by motorized pump permitting the sediment at the bottom to be easily removed.
- The well is deepened until clean sand is reached, then a layer of coarse gravel is put on the bottom.
- A strong solution of chlorine (e.g. 4 L of 5% available chlorine solution) is then put into the water and the well is allowed to recharge. It may take 2-3 days for the chlorine to dissipate depending upon the volume of water treated and the amount of organic matter that still remains in the well.

Notes: Water taken from an open well, even though it may look clean, should be considered contaminated and should be boiled or disinfected. This is due to the fact that you cannot prevent contaminants from getting into the well. It is impossible to effectively disinfect an open well that has not been lined. The interior surface of the well must be smooth and easily cleanable.

Chlorination of a well can be done by the tube method. This technique allows a slow release of chlorine over a longer period of time than the dosing method. Such devices can leave a chlorine residual for sometimes 2-3 weeks depending upon the volume of the well, its use and the amount of free chlorine available. Water and health authorities should be consulted for the proper design and maintenance procedures of these types of chlorinators.

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5.3 Storage

Very often we forget that drinking water can become contaminated not only at the source or in the piped distribution system, but also during storage at home. If we want to prevent the transmission of waterborne diseases, then we must safeguard against the re-contamination of water in the home.

A two-component strategy to prevent diarrhoeal diseases in households that lack a safe continuous supply of water should consist of:

- Water disinfection immediately after collection (or filtration) ; and
- Safe storage: storage in vessels designed to prevent re-contamination.

Successful implementation of this strategy will require an educational campaign stressing the role of contaminated water and domestic hygiene in disease transmission.

Numerous studies around the world have found higher concentrations of disease - causing organisms in water storage containers in the home than in the water sources. During the cholera epidemic in Peru, for example, concentrations of faecal coliforms in several peri-urban districts were found to be up to a thousand times higher in storage vessels than in municipal water taps. Similarly, studies in Asia and the Middle East have detected the presence of organisms like *V. cholerae*, enterotoxigenic *E. coli* and *Ascaris* in household water containers while these organisms were absent in water samples from the taps and wells where the water was originally collected from.

Re-contamination of water in storage containers usually comes from hands and/or objects (eg. Cups) that are introduced to get water out. Narrow opening preventing the introduction of objects reduces significantly the risk of re-contamination. Covering the opening of the container is another important safeguard. The material of the storage vessel itself has also been found to influence the survival of bacteria in the stored water.

For example, in inoculation experiments with African storage containers, *V. cholerae* survived as long as:

- 7 days in clay pots,
- 22 days in plastic containers, and
- 27 in metal drums.

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

Not surprisingly, survival times did not exceed one day with chlorine levels of 0.2 mg/l. or greater in any type of container.

For safe storage, a water storage vessel should have the following qualities:

- Be constructed of a material that is durable, light weight, non-oxidizing, and easy to clean.
- Hold an appropriate standard volume (e.g. 20 L) for ease of dosing disinfectant.
- Have a stable base and a sturdy, comfortable handle for easy carriage.
- Have a single opening 5 to 8 cm in diameter with a strong, tight fitting cover-this dimension makes it easy to fill container and add disinfectant but difficult to insert hands or utensils.
- Have a non-rusting, durable, cleanable spigot for extracting water.
- Allow air to enter as water is extracted.
- Have volume indicators (for disinfectant dose calculations).

Traditional water storage vessels that have been modified to reduce the risk of contamination: (A): Egyptian zir; (B) Vegetable oil container from Zambia; (C) Cantero from El Salvador; (D) Sorai from India; (E) Tin bucket from Malawi; (F) Plastic container designed by the US Centre for Disease Control and Prevention.

Illustrated directions on how to fill the water storage vessel, add proper amount of disinfectant, and remove water for drinking.

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6. Training Evaluation & Closing

Evaluation Test for Trainees:

An oral test should be given at the end. The trainees are told days in advance that they will have an oral test. Depending on the group, they are split into 4-5 teams (e.g. ward against ward). Modest prizes can be awarded to the winning team (the one that answers most questions correctly). These prizes can be for example, a pen, chocolates, or anything that is inexpensive and valued by the participants. Address questions to the different groups in order (first question to group A, second question to group B, and so on) BUT give only 30 seconds for each group to give the correct answer. If no correct answer by group A is given, then go to B, then to C, then D. If none gives the correct answer, open the answer to everyone. Usually someone who was not in the group decision will respond. Give the point to the person's team. For some questions, there may be 4-5 correct answers, still give 30 seconds and the next group may give you the rest, or the next group. Always score one point for correct answer. Soon you will find that the groups not being asked will try to work out the answer because they know they may have a chance to answer the question. This exercise can provide an excellent review and everyone listens and learns. About 30 questions should be made, totaling 40-50 points.

Feedback:

Finally, it is always good to ask the participants for feedback on the training: what they liked best and least; what they had the most difficult with and why, etc.....

Evaluating the WQM Program:

This particular training program is meant to be practical in nature. It is a "hands-on training" program which not only explains what to do, but allows the trainees to do it themselves. Therefore, to evaluate the effectiveness of the training it is necessary to assess skills. This may be done during regular visits by trainers by comparing the records with what actually exists and observing what is being done. Sampling of sites should also be made during these visits. It is recommended that for about 10% of the samples, replicates with Coliplates or Colistrips be carried out as quality control, as well as to compare with the performance of the H₂S test.

If routine maintenance or repair work or disinfection is not being done, then the project team is questioned to determine the reasons-which may be technical, administrative or financial. If however it is due to lack of knowledge, then follow-up training must be done.

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At this time a decision can be made to either explain a procedure or (if complicated) to wait until an opportune time to train several WQM technicians at the same time.

Having wards assist other wards in the training program is another way to verify if the information presented during the training was understood. It is also a good way to motivate new communities to become involved as they will be proud to show what they have done.

Another type of evaluation that might be done is a health impact study. If the water quality monitoring and control program is functioning as intended, then diarrhoeal disease rates due to water borne diseases should decline over time. This type of evaluation however, would require the services of an epidemiologist to design the research protocol.

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WATER QUALITY MONITORING TRAINING PROGRAM

ANNEXES:

Trainee Handouts

Community Mapping

Water Use Survey

Transact Walk

The H₂S Test

Coliplate™ & Colistrip™ Test Guide

[illegible]

ANNEX: Water Use Survey (Matrix voting example)

SOURCE					
USE					
WASH BODY					
WASH CLOTHES					
GARDEN					

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ANNEX: Transect Walk

Sample Recording Form: Water Source/System - Reliability - access - usage - water quality

Ward: _____

Date: _____

Water Source/System Location: _____

Recorded by: _____

Item	Yes	No	Don't know	Disagreement in group	Remarks
Does water source runs dry/low during the year? (Circle condition) If yes, when and how long: _____					
Are households in the surroundings using this water source? If yes, what are the uses? (Circle uses) - Drinking? - bathing? - washing? - cooking? - rituals? - animals? - gardening/irrigation? Other: _____?					
Is it NOT used for some purposes? Which: _____ Why: _____					
Is water source accessible to all household in neighborhood? If not, why? _____					
How many households have no access?					
Is the quality of the water acceptable? If no, what is wrong with it? _____					
Is the water quality tested regularly? If yes, for what? _____					
Are results of testing reported to the community?					
Is action taken when results are poor? What: _____					
Is there someone responsible for the upkeep of the water source? Has she/he made any repairs? If yes, what?: _____					
Has this water point ever broken down? If yes, what was the quickest repair made? What was the longest one?					

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ANNEX: Transect Walk Sample Scoring & Recording Form (Sanitation)

Ward :

Date:

Location:

Recorded by:

Score: 2=good 1=average 0=poor	Scores for Ward Water Points: (P): Public (H): Household									
	Health Clinic	# ()	# ()	# ()	# ()	# ()	# ()	# ()	# ()	# ()
1. Latrine/Toilet in working order										
2. Latrine/Toilet used regularly										
3. Pit/facility built properly										
4. Pit/facility safely located*										
5. Superstructure offers privacy of use										
6. Cover present & over hole/water in water seal										
7. No excreta visible on floor/walls/in pan										
8. Water and soap/soap substitute at or near facility										
9. No stagnant water in inside floor nor outside facility										
10. No adult/child excreta around house or yard										
Total score (out of 10)										

1. Latrine pit or septic tank at least 15 meters away and downstream from water source.
2. If no latrine or toilet, or facility never used, score a "0" in total score for that household.

NB. Additional sheets can be used to rank more latrines/toilets depending on project & village size and context.

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ANNEX: Transect Walk Sample Scoring & Recording: Reliability & Quality of Water Sources

Ward :

Date:

Location:

Recorded by:

Score: 2=good 1=average, 0=poor	Scores for Ward Water Points: (P): Public(H): Household									
Reliability (R)	Health Clinic	# ()	# ()	# ()	# ()	# ()	# ()	# ()	# ()	# ()
Water Point Source Reliability (Overall) RO										
Reliability (Women's needs) RW										
Reliability (Men's needs) RM										
Total (R) (Out of 6)										
Water Quality Setting (WQS)										
Water Quality according to users: WQU										
Water Point Source Protection: WP										
Water Testing: WT										
Total (WQS) (Out of 6)										

Reliability (R)	RO: Source (Overall)	RW: Women's needs	RM: Men's needs
2 = good	ALWAYS ENOUGH water to meet basic needs	CAN ALWAYS fulfill all water needs from source	CAN ALWAYS fulfill all water needs from source
1 = average	Source never goes dry but NOT ENOUGH water for SEVERAL DAYS yr	NOT ENOUGH water to meet needs for SEVERAL DAYS during year	NOT ENOUGH water to meet needs for SEVERAL DAYS during year
0 = poor	EXTENDED PERIODS of time with NO WATER available at the source	NOT ENOUGH water to meet needs for SEVERAL WEEKS during year	NOT ENOUGH water to meet needs for SEVERAL WEEKS during year
Water Quality Setting (WQS)	WQU: WQ according to users	WP: Source Protection	WT: Water Testing
2 = good	Always of good quality	Source is protected and water is treated	Regular testing for bacteria people informed of results actions taken when contamination is detected
1 = average	Poor quality only during certain times of year (list problems & times):	Source is protected or water is treated	Water sometimes tested for bacterial but results not reported to people
0 = poor	Usually of poor quality (list problems):	Sources with no protection of protection treatment not working	Water testing never done

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ANNEX: The H₂S Test

Hygiene and Safety Precautions:

The following precautions are essential, not only to avoid contamination of the samples and culture medium, but also to avoid the danger of infection to staff conducting the microbiological tests.

- Keep finger nails short: use hair and beard protection if necessary.
- Wash hands with soap and warm water before and after microbiological examinations as well as after use of the toilet.
- Do not touch your mouth or eyes while working with the test bottles.
- If you cough or sneeze, cover your mouth. Clean and dry your hands before continuing work with the microbiological preparations or examinations.
- Keep your work clothes clean.
- If you feel sick you should not work with the tests.
- Do not eat, drink or smoke in the working area.

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- Washing the surface with detergent and water removes dust but does not eliminate all germs. It is essential to disinfect the working area. A cheap solution to disinfect all surfaces after they have been washed is to use bleach or chlorine (home disinfectant with 5.25% of available chlorine) mixed with water.

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- The chlorine solution must be left to rest for at least 30 seconds on the surface to be disinfected to ensure this is well done.

In case of accidents or spills:

- Immediately clean all dangerous debris.
- Place all broken glassware in a closed garbage container with a plastic bag lining.
- If a contaminated (or potentially contaminated) sample is spilled or its container is broken, disinfect the area with a chloride solution.
- Use rubber gloves during cleaning procedures.
- It is very important not to touch your eyes or mouth with our hands.
- If you get hurt, go see a nurse or a physician.

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H₂S media:


Bacteriological peptone.....	40.0 g
Dipotassium hydrogen phosphate	3.0 g
Ferric ammonium citrate.....	1.50 g
Sodium thiosulphate	2.00 g
Teepol	2.0 ml.
Water, distilled or boiled tap.....	100.0 ml.

H₂S media volumes:


The H₂S test allows to estimate different degrees of contamination by using different volumes of water samples with higher or lower amounts of H₂S media impregnating the paper strips. The following are recommended variations.

Test Sample Volume	Volume of H ₂ S media added to paper strip	Recommended Use
1 ml.	0.5 mL	untreated waters or waters with suspected contamination
10 ml.	0.5 mL	untreated waters
20 ml.	1 ml.	treated waters
100 ml.	2.5 mL	treated waters in well maintained systems


Marking volumes in test bottles & tubes:

- 1. Prepare “marking guide tubes & bottles” with graduated cylinder or pipettes
- 2. Use guides to mark test containers before use


1 mL tests: tube with 9 & 10 mL markings
10 mL tests: tube with 10 mL markings


20 mL tests: bottle with 20 mL markings
100 mL tests: bottle with 100 mL markings

Inexpensive field incubator using insulated ice-chest


Figure

H₂S Field Monitoring Qualitative Procedure

1 mL positive for H₂S = 100 + coliforms / 100mL
 1 mL negative for H₂S = < 100 coliforms / 100mL

10 mL positive for H₂S = 10 + coliforms / 100 mL
 10 mL negative for H₂S = < 10 coliforms / 100 mL

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A "control" (10 mL boiled water tube) must be used for each new source of boiled water used. If the "control" gives a positive result (black coloration), this indicates that the 1 mL samples to which this boiled water was added can not be considered valid. It may also indicate an inadequate procedure resulting in the contamination of the test tubes. In this case carefully repeat the water sampling and its examination.

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Community : _____

[illegible]

Data Collections sheets - H₂S Test (10 and 1 mL sample volumes)

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Data Collections Sheets - H₂S Test (Traditional Method)

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ANNEX: Coliplate™ & Colistrip™ Test Guide

Coliplate™ MPN Procedures for Total Coliforms and *E.coli*

The Coliplate™ is a convenient test for the quantitative measure of total Coliforms and *E.coli*. The test is designed to meet regulatory guidelines for surface water and raw drinking water, as well as recreational waters and wastewater. The full Coliplate enables quantification of Coliforms and *E.coli* densities ranging from <3 to 2400 colony forming units (cfu)/100 mL in a single test without dilution. The distinctive blue/green coloration of positive (Coliforms) tested samples enables analysis of brownish, turbid or rust containing waters.

The Coliplate test utilizes the proven X-Gal and MUGi techniques to detect viable Coliforms and/or *E.coli*. The Coliplate contains selective media to provide nutrients to stimulate the growth of Coliforms and *E.coli*. The media also contains inducers and chromogenic/fluorogenic substrates. These substrates react with enzymes specific to coliform (beta-D-galactosidase) and *E.coli* (beta-D-glucuronidase) to provide colour change to blue/green (Coliforms) and fluorescence, (*E.coli*).

Test results are recorded after 20-24 hours of incubation with the appearance of a blue/green colour in the wells containing Coliforms. *E.coli* can be detected by the appearance of fluorescence under a long wavelength UV light, in some or all the wells a blue/green colour. The Coliplate can be used to confirm H₂S paper strip results and establish *E.coli* densities.

Coliplate™ - Whole Plate Procedure:

1. Using a pipette fill up all 96 wells of the microplate; or dip the whole Coliplate into the sample and fill up all 96 wells. This last procedure works with well samples, where the water is collected in a bucket.
2. Cover microplate with lid; insert microplate into zip-loc bag (or other plastic bag which can be sealed) and incubate at 35° C to 20-24 hours.
3. After incubation, count the number of wells that have turned blue/green. Refer to 96 well MPN Table to determine the Most Probable Number for total Coliforms in 100 mL sample.
4. Place long wave UV (366 nm) light over the microplate and count the number of wells that have a blue fluorescence. Refer to 96 well MPN Table to determine the Most Probable Number for *E.coli* in 100 mL sample.
5. For sample which are very polluted (i.e. 1 mL H₂S test is positive in 24 hours or less) it is recommended that only 1/2 of the microplate is used (48 wells). Thus 2 samples can be tested on the same plate. For MPN Total Coliforms and *E.coli* refer to the 48 well MPN.

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
Coliplate™ - Two Strip Quick Screening Procedure:

1. Dip the double strip unit into the water sample until all wells are filled up completely.
2. Pull the strip unit back gently and place inside a plastic bag and seal. or
3. Using a pipette fill up 2 rows of wells (16 wells) of a 96 well microplate. Label. with this procedure 6 water samples can be tested with on 96 well microplate.
4. Incubate at 35^o C for 20-24 hours.
5. Count the number of wells that turned blue/green. Refer to 16 well Microplate MPN Table to determine the Most probable Number for Total Coliforms in 100 mL sample.
6. Observe under long wave UV (366 nm) light and count the number of wells that have a blue fluorescence. Refer to 16 well Microplate™ MPN Table for *E.coli* in 100 mL sample.

No. Of wells giving positive reaction	MPN per 10 mL sample	No. of wells giving positive reaction	MPN per 10 mL sample	No. of wells giving positive reaction	MPN per 10 mL sample	No. of wells giving positive reaction	MPN per 10 mL sample
0	<16						
1	16	5	94	9	213	13	469
2	33	6	119	10	255	14	619
3	52	7	146	11	307	15	938
4	72	8	177	12	375	16	>938

Field Monitoring Quantitative Coliplate™-400 Procedure for E. Coli and Total Coliforms

Blue colour = coliform
fluorescing blue colour = *E. Coli*

Full plate range
< 3 - > 2424
with 0.40 mL addition
per cup

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