


A Whole Cell Biosensor (*GlnLux*) to Measure Symbiotic Nitrogen Fixation (SNF) in Legumes

Abstract

Measuring symbiotic nitrogen fixation (SNF) of large numbers of legume plants is a challenge in terms of time, labor, and cost. Here we present a whole cell biosensor to measure SNF efficiently in legumes. A whole cell biosensor (*GlnLux*) for glutamine (Gln) was constructed by transforming a bacterial Gln auxotroph with a constitutive lux reporter. Extracts from single legume leaf punches were sufficient for measuring SNF by the *GlnLux* method, wherein leaf extracts are incubated with *GlnLux* and luminescence is measured using a luminometer in 96 well plates. This method has potential to measure SNF in amide exporting legumes (lentil) as well as ureide exporting legumes (cowpea). Further the *GlnLux* method can be used to identify the effect of different rhizobia and plant genotypes on SNF. Interestingly the *GlnLux* method permitted localization of active sites of SNF in nodulated legume roots. For imaging of Gln in whole plant organs, tissues were freeze-thawed to cause Gln leakage, placed on agar pre-embedded with *GlnLux* (*GlnLux* agar), and then imaged using a photon capture camera.

Construction of *GlnLux* Whole Cell Biosensor



GlnLux Whole Cell Biosensor Assays to Measure SNF




Figure 1 | (a) **Overview of concept:** Rhizobia inside the nodules (bacteroids) convert atmospheric nitrogen into NH_4^+ , which is assimilated into Gln. In *in vitro* bioassays, Gln is taken up by *GlnLux* biosensor cells, allowing them to multiply, replicate the lux plasmid and release photons. (b) **96-well liquid assay:** Leaf punches are sampled from the youngest full expanded leaves and immediately frozen in liquid nitrogen. Individual frozen leaf punches are ground using a micropestle with silica sand in a microcentrifuge tube placed on ice with protease inhibitor cocktail (PIC). Plant extracts are centrifuged, and the supernatant is transferred to a microcentrifuge tube on ice. White opaque 96-well reader plates are loaded with M9 minimal medium to which each plant extract is added. Finally, an aliquot of 14 h-Gln-depleted *GlnLux* culture is added to each well. Plates are sealed with sterile film, mixed, and incubated at 37 °C without shaking. For lux quantification, plates are read in a luminometer. (c) **Solid assay (tissue imaging):** Legume roots are frozen in -80 °C and thawed at room temperature for 1 min. Roots are placed in contact with agar containing *GlnLux* cells (*GlnLux* agar) and tissues are pressed down. Inverted plates are imaged at the zero time point, then incubated at 37 °C with hourly imaging with 300-600 s exposure times using a ChemiProHT luminescence imaging system.

GlnLux Distinguishes Amide-Exporting Plants Inoculated with Diverse Rhizobia

GlnLux Liquid Assay with leaf punches




Figure 2 | (a) Lentil plants inoculated with four different strains of *R. leguminosarum* bv. *viciae* (wild type: 3841, VF39, 248; mutant 17-B) and uninoculated control. (b) Mutant strain (17-B) produced smaller nodules compared to the plants inoculated with wild type strains. (c) Shoot Gln content of lentil plants tested using the *GlnLux* assay (RLU, relative light unit). (d) Shoot total N% of lentil plants. (e) The % of nitrogen derived from atmosphere (%Ndfa) of lentil plants calculated using ^{15}N isotope dilution technique. (f) The correlation between *GlnLux* and %Ndfa.

GlnLux Solid Assay with Roots




Figure 3 | Luminescence imaging of pure Gln liquid standards and lentil root Gln using *GlnLux* agar. (a-c) Agar discs were prepared with different concentrations of Gln (0, 10^{-7} - 10^{-2} M) and placed on *GlnLux* agar. The opposite agar surface was then imaged using a photon capture CCD camera. (d-f) Four-week-old lentil roots inoculated with four different strains of *R. leguminosarum* bv. *viciae* (VF39, 248, 3841, 17-B) or control (un-inoculated) were freeze-thawed to cause Gln leakage, and then placed on *GlnLux* agar. The opposite agar surface was then imaged using a photon capture CCD camera. (a,d) The light image, (b,e) white lux image, and (c,f) false-colored lux images after 3 hours of incubation and 300 s exposure.

GlnLux Distinguishes Ureide-Exporting Plants Inoculated with Diverse Rhizobia




Figure 4 | (a) Cowpea plants inoculated with *Bradyrhizobium yuanmingense* (TTC9, TSC10), *B. japonicum* (DTB4), *B. elkanii* (DTC9), or an un-inoculated control. (b) Number of nodules per plant. (c) Shoot Gln content of lentil plants tested using *GlnLux* leaf punch assay (RLU, relative light unit). (d) Shoot total N% of lentil plants. (e) The %N derived from the atmosphere (%Ndfa) in lentil plants calculated using ^{15}N isotope dilution. (f) The correlation between *GlnLux* and %Ndfa.

Active Sites of Nitrogen Fixation can be Visualized Using *GlnLux*




Figure 5 | Split root experiment with soybean plant. (a-b) One side was inoculated with *B. japonicum* 110 (WT), while the other side was inoculated with 510 (MT). (c-e) Luminescence *in vivo* imaging of Gln in soybean roots using *GlnLux* agar. Freeze-thawed roots were placed on *GlnLux* agar. The opposite agar surface was then imaged using a photon capture CCD camera. (c) The light image, (d) white lux image, and (e) false-colored lux images after 2 hours of incubation and 600 s exposure.

Conclusions

- GlnLux* biosensor is a new method to measure SNF in both ureide- and amide-exporting legumes.
- GlnLux* 96-well liquid assay uses single leaf punches to measure relative SNF in plants growing without exogenous N, making it a rapid, low-cost, high throughput screening method.
- GlnLux* solid assay permits visualization of active sites of nitrogen fixation.
- GlnLux* method can be used to screen plants inoculated with different rhizobia strains to optimize SNF.
- GlnLux* method is being tested for its ability to distinguish plant genotypes that vary in SNF.

References

Tessaro, M.J., Soliman, S.S.M. and Raizada, M.N. 2012. Bacterial Whole-Cell Biosensor for Glutamine with Applications for Quantifying and Visualizing Glutamine in Plants. *Appl. Environ. Microbiol.* **78**, 604–606.
Thilakarathna, M.S. and Raizada, M.N. 2015. Bacterial Whole-cell biosensor (*GlnLux*) to Measure Symbiotic Nitrogen Fixation in Legumes. (in preparation)

Acknowledgements

We thank Jaelyn Clark and Caleb Niemeyer (University of Guelph) for their significant help in data collection. We acknowledge Dr. Chris Yost and Dr. Dinah Tambalo (University of Regina), Prof. Yamakawa (Kyushu University) for providing rhizobia strains, and Kamal Khadka (University of Guelph) for providing seeds. Special thanks to Lisa Smith for preparing graphic illustrations.